
FINESSE[®] *Ultra*

Introducer Stylet

Instructions for Use

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BIOPSY SYSTEMS

INSTRUCTIONS FOR USE

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

Refer to the **FINESSE® ULTRA Biopsy Driver and Probe Instructions for Use** for additional important information.

Single-use only.

A. Device Description

The **FINESSE® ULTRA Introducer Stylet** is part of the **FINESSE® ULTRA Breast Biopsy System**. It is composed of a stylet, hub, and protective sheath and is only to be used with the **FINESSE® ULTRA Breast Biopsy System**. The components of the system are designed to operate safely when used together for diagnostic sampling of tissue during a biopsy procedure.

B. Indications for Use

The **FINESSE® ULTRA Breast Biopsy System** is indicated to obtain tissue samples from the breast or axillary lymph nodes for diagnostic analysis of breast abnormalities. The instrument is intended to provide breast tissue for histologic examination with partial or complete removal of the imaged abnormality. The extent of histologic abnormality cannot be reliably determined from its mammographic appearance. Therefore, the extent of removal of the imaged evidence of an abnormality does not predict the extent of removal of a histologic abnormality e.g., malignancy. When the sampled abnormality is not histologically benign, it is essential that the tissue margins be examined for completeness of removal using standard surgical procedures.

C. Contraindications

1. The **FINESSE® ULTRA Breast Biopsy System** is for diagnostic use only, NOT for therapeutic use.
2. The **FINESSE® ULTRA Breast Biopsy System** is contraindicated for those patients where, in the physician's judgment, there is an increased risk of complications associated with percutaneous removal of tissue samples.

D. Warnings

1. **Patients who may have a bleeding disorder or who are receiving anticoagulant therapy may be at increased risk of complications.**
2. **As with any biopsy instrument, there is a potential for infection.**
3. **Do not reuse the FINESSE® ULTRA Introducer Stylet. Reusing the stylet bears the risk of cross-patient contamination as the residue of biological material can promote the contamination of the stylet with pyrogens or microorganisms which may lead to infectious complications.**
4. **Do not resterilize the FINESSE® ULTRA Introducer Stylet. After resterilization, the sterility of the stylet is not guaranteed because of an indeterminable degree of potential pyrogenic or microbial contamination which may lead to infectious complications.**
5. **The FINESSE® ULTRA Introducer Stylet must only be used with the FINESSE® ULTRA Breast Biopsy System.**
6. **All breast biopsies should be performed under imaging guidance to confirm the stylet position relative to the target region to be sampled and to help mitigate the occurrence of a false negative biopsy. The FINESSE® ULTRA Introducer Stylet is intended for use with ultrasound imaging only.**

E. Precautions

1. Do not use the **FINESSE® ULTRA Introducer Stylet** without the **FINESSE® ULTRA Biopsy Probe** integrated coaxial cannula.
2. The **FINESSE® ULTRA Breast Biopsy System** should only be used by a physician trained in its indicated use, limitations, and possible complications of percutaneous needle techniques.

F. Potential Complications

1. Potential complications are those associated with any percutaneous removal/biopsy technique for tissue collection. Potential complications are limited to the region surrounding the biopsy site and include hematoma, hemorrhage, infection, a non-healing wound, pain and tissue adherence to the biopsy probe while removing it from the breast.
2. As per routine biopsy procedures, it may be necessary to cut tissue adhering to the introducer stylet while removing it from the breast.

G. Equipment Required

The following equipment is required for a biopsy procedure:

- Appropriate imaging modality and accessories
- **FINESSE® ULTRA Biopsy Driver**
- **FINESSE® ULTRA Biopsy Probe**
- **FINESSE® ULTRA Introducer Stylet**
- Surgical gloves and drapes
- Local anesthetic
- Scalpel
- Other equipment as necessary

H. Directions for Use

As a reference, please refer to the directions for use for both the probe and driver.



1. Load the probe into the driver in accordance with the probe and driver IFU. Leave the protective sleeve on the probe needle. Ensure that the probe remains sterile.

2. Remove the integrated coaxial cannula from the probe. (Figure 1)

3. Remove the introducer stylet from the package and discard the needle protector. Ensure the introducer stylet remains sterile.


4. Install the introducer stylet into the coaxial cannula. (Figure 2) Ensure the coaxial cannula is locked onto the introducer stylet hub. (Figure 3)

5. Prepare the biopsy site using appropriate aseptic techniques and adequate local anesthesia. Appropriate precautions should be followed for any percutaneous procedure (personal protective gear such as gloves and glasses/shields, etc. should be worn)

6. Make a small incision or puncture the skin with a scalpel.

7. The introducer stylet features raised bumps to indicate the correct way up. Grip the introducer stylet such that the thumb rests on the raised bump surface. (Figure 4)

8. Using ultrasound guidance, insert the introducer stylet through the incision and position the tip to the appropriate depth.
9. With the raised bump surface oriented facing up, remove the stylet from the integrated coaxial cannula. The integrated coaxial cannula will retain a track to the biopsy site.
10. Insert the probe into the integrated coaxial cannula and ensure it is locked into place.
11. Perform the tissue sampling procedure as described in the probe instructions for use.

NOTE: The introducer stylet, probe, coaxial cannula, and sample container are for single use only. After use, these products may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

I. How Supplied

The FINESSE® ULTRA Introducer Stylet is supplied sterile for single patient use only.

J. Warranty

Bard Peripheral Vascular warrants to the first purchaser of this product that this product will be free from defects in materials and workmanship for a period of one year from the date of first purchase and liability under this limited product warranty will be limited to repair or replacement of the defective product, in Bard Peripheral Vascular's sole discretion or refunding your net price paid. Wear and tear from normal use or defects resulting from misuse of this product are not covered by this limited warranty.

TO THE EXTENT ALLOWABLE BY APPLICABLE LAW, THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL BARD PERIPHERAL VASCULAR BE LIABLE TO YOU FOR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM YOUR HANDLING OR USE OF THIS PRODUCT.

Some states/countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your state/country.

An issue or revision date and a revision number for these instructions are included for the user's information on the last page of this booklet. In the event 36 months have elapsed between this date and product use, the user should contact Bard Peripheral Vascular, Inc. to see if additional product information is available.

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Copyright ©2011 C. R. Bard, Inc. All rights reserved.
Copyright ©2011 C. R. Bard, Inc. Tous droits réservés.
Copyright ©2011 C. R. Bard, Inc. Alle Rechte vorbehalten.
Copyright ©2011 C. R. Bard, Inc. Tutti i diritti riservati.
Copyright ©2011 C. R. Bard, Inc. Todos los derechos reservados.
Copyright ©2011 C.R. Bard, Inc. Alle rechten voorbehouden.