
ULTRACLIP® II Tissue Marker

Instructions for Use

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BIOPSY SYSTEMS


Instructions for Use

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

A. General Information and Device Description:

The ULTRACLIP® II tissue marker is a sterile, single use device comprised of a disposable introducer, and a metal implantable tissue marker clip. The introducer needle has 1 cm reference marks and an ultrasound enhancement on the distal end to aid in needle placement. A safety switch on the introducer helps prevent premature deployment of the marker (See Figure 1).

Figure 1

1. 1 cm reference marks
2. Introducer needle
3. Safety switch
4. Plunger
5. ULTRACLIP® II tissue marker

The marker, located inside the distal end of the introducer needle, is made of Titanium, Inconel® 625 or BioDur® 108.

Product Code	Gauge & Length	Marker Material	Marker Shape
861017	17g x 10cm	Titanium	Ribbon 
861217	17g x 12cm	Titanium	Ribbon 
863017	17g x 10cm	Titanium with PVA	Ribbon 
865017	17g x 10cm	Titanium	Ribbon 
865517	17g x 15cm	Titanium	Ribbon 
862017	17g x 10cm	Inconel® 625	Wing 
864017	17g x 10cm	BioDur® 108	Coil 

Product Code #863017 has a PVA polymer pellet interwoven with the titanium clip to enhance viewing of the clip under ultrasound. Please note this polymer is non-absorbable.

Product codes #865017 & #865517 are MRI Conditional and may be placed under MRI guidance.

B. How Supplied:

The ULTRACLIP® II device is supplied sterile and non-pyrogenic unless the package has been opened or damaged.

Sterilized using Ethylene Oxide. For single use only. Do not reuse. Do not resterilize.

C. Indications for Use:

The ULTRACLIP® II tissue marker is intended for use to attach to soft breast tissue at the surgical site during an open surgical breast biopsy or a percutaneous breast biopsy to radiographically mark the location of the biopsy procedure.

D. Contraindications for Use:

This device is not intended for use except as indicated above.

E. Warnings:

1. The ULTRACLIP® II tissue marker is not recommended for use in patients with breast implants.
2. As with any foreign object implanted into the body, potential adverse reactions are possible. It is the responsibility of the physician to evaluate the risk/benefit prior to the use of this device.
3. This device has been designed for single use only. Reusing this medical device bears the risk of cross-patient contamination as medical devices – particularly those with long and small lumina, joints, and/or crevices between components – are difficult or impossible to clean once body fluids or tissues with potential pyrogenic or microbial contamination have had contact with the medical device for an indeterminate period of time. The residue of biological material can promote the contamination of the device with pyrogens or microorganisms which may lead to infectious complications.
4. Do not resterilize. After resterilization, the sterility of the product is not guaranteed because of an indeterminable degree of potential pyrogenic or microbial contamination which may lead to infectious complications. Cleaning, reprocessing and/or resterilization of the present medical device increases the probability that the device will malfunction due to potential adverse effects on components that are influenced by thermal and/or mechanical changes.

F. Precautions:

1. The marker must only be deployed by depressing the plunger at the proximal end of the device. Do not attempt to deploy the marker by pressing or pulling on the safety switch.
2. This product should only be used by a physician who is completely familiar with the indications, contraindications, limitations, typical findings and possible side effects of tissue marker placement.
3. After use, this product may be a potential biohazard. Handle and dispose of in accordance with acceptable medical practice and applicable local, state and federal laws and regulations.

G. Potential Complications:

Potential complications of marker placement may consist of hematoma, hemorrhage, infection, adjacent tissue injury and pain.

H. Equipment Required:

The following equipment is required to place the marker:

- Appropriate imaging modality and accessories
- Surgical gloves and drapes
- Local anesthetic
- Other equipment as necessary

I. Directions for Use:

1. Prepare the site as required. Adequate anesthesia should be administered, as required.
2. Inspect the package and product for damage. If undamaged, open the package and transfer the ULTRACLIP® II device onto the sterile field utilizing aseptic technique.

3. Locate the target area for deployment using the appropriate imaging technique.
4. Insert the introducer needle into the breast, directing it to the target. Use the 1 cm reference markings to position the needle point just proximal to the target. **NOTE:** The BARD® logo and the safety switch line up with the level of the needle to help in needle placement.
5. Confirm needle placement with the appropriate imaging technique. If necessary, reposition the needle and reconfirm placement.
6. Release the safety switch by moving the switch to line up with the longitudinal slot on the device (**Figure 2-a**) or by rotating the plunger at the proximal end of the device (**Figure 2-b**).

Figure 2-a




Figure 2-b

7. Firmly depress the plunger with the thumb or index finger until a firm click is heard or felt (**Figure 3**). (Do not attempt to deploy the marker by pressing or pulling on the safety switch.)

Figure 3

8. Remove the introducer and confirm placement of the marker using the appropriate imaging technique.

**MR Conditional Information:**

Non-clinical testing has demonstrated that the ULTRACLIP® II tissue marker is MR Conditional. It can be scanned safely under the following conditions:

- Static magnetic field of 3-Tesla or less
- Spatial gradient field of 720-Gauss/cm or less
- Maximum MR system reported whole-body-averaged specific absorption rate (SAR) of 3-W/kg for 15 minutes of scanning.

In non-clinical testing, the marker produced a temperature rise of less than +0.4°C at a maximum MR system reported whole body averaged specific absorption rate (SAR) of 3-W/kg of 15 minutes of MR scanning in a 3-Tesla, Signa MR system (Excite platform; G3.0-052B software, GE Healthcare, Milwaukee, WI).

Product codes **865017** and **865517** can be placed and imaged under MR. In non-clinical testing, the loaded introducer of these product codes produced a temperature rise of less than +1.3°C at a maximum MR system reported whole body averaged specific absorption rate (SAR) of 3-W/kg of 15 minutes of MR scanning in a 3-Tesla, Signa MR system (Excite platform; G3.0-052B software, GE Healthcare, Milwaukee, WI).

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the ULTRACLIP® II tissue marker. Therefore, it may be necessary to optimize MR imaging parameters for the presence of this metallic implant.

Warranty

Bard Peripheral Vascular warrants to the first purchaser of this product that this product will be free from defects in materials and workmanship for a period of one year from the date of first purchase and liability under this limited product warranty will be limited to repair or replacement of the defective product, in Bard Peripheral Vascular's sole discretion or refunding your net price paid. Wear and tear from normal use or defects resulting from misuse of this product are not covered by this limited warranty.

TO THE EXTENT ALLOWABLE BY APPLICABLE LAW, THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL BARD PERIPHERAL VASCULAR BE LIABLE TO YOU FOR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM YOUR HANDLING OR USE OF THIS PRODUCT.

Some states/countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your state/country.

An issue or revision date and a revision number for these instructions are included for the user's information on the last page of this booklet. In the event 36 months have elapsed between this date and product use, the user should contact Bard Peripheral Vascular, Inc. to see if additional product information is available.

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BARD and ULTRACLIP are trademarks and/or registered trademarks of C. R. Bard, Inc. or an affiliate. All other trademarks are the property of their respective owners.
Copyright ©2013 C. R. Bard, Inc. All rights reserved.
Copyright ©2013 C. R. Bard, Inc. Tous droits réservés.
Copyright ©2013 C. R. Bard, Inc. Alle Rechte vorbehalten.
Copyright ©2013 C. R. Bard, Inc. Tutti i diritti riservati.
Copyright ©2013, C. R. Bard, Inc. Reservados todos los derechos.