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
BIOPSY SYSTEMS

Instructions for Use

Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

A. General Information and Device Description:

The BARD® DUALOK® Breast Lesion Localization Wire is a two-component device consisting of a needle cannula with integral female luer lock hub and a localization wire with dual locking elements at the distal tip. The BARD® DUALOK® breast lesion localization wire is a disposable single patient use device and is available in various lengths.

B. How Supplied:

The product is supplied sterile and nonpyrogenic unless the package has been opened or damaged. **Sterilized using Ethylene Oxide. For single use only.**

Do not reuse. Do not resterilize.

C. Indications for Use:

The BARD® DUALOK® Breast Lesion Localization Wire is intended for use as a guide for directing surgical biopsies.

D. Contraindications:

None known.

E. Warnings:

1. The BARD® DUALOK® Breast Lesion Localization Wire has been designed for single use only. Reusing this medical device bears the risk of cross-patient contamination as medical devices – particularly those with long and small lumina, joints, and/or crevices between components – are difficult or impossible to clean once body fluids or tissues with potential pyrogenic or microbial contamination have had contact with the medical device for an indeterminable period of time. The residue of biological material can promote the contamination of the device with pyrogens or microorganisms which may lead to infectious complications.
2. Do not resterilize the BARD® DUALOK® Breast Lesion Localization Wire. After resterilization, the sterility of the product is not guaranteed because of an indeterminable degree of potential pyrogenic or microbial contamination which may lead to infectious complications. Cleaning, reprocessing and/or resterilization of the present medical device increases the probability that the device will malfunction due to potential adverse effects on components that are influenced by thermal and/or mechanical changes.

NOTE: After use, this product may be a potential biohazard. Handle and dispose of in accordance with acceptable medical practice and applicable local, state, and federal laws and regulations.

F. Precautions:

1. This device should only be used by a physician trained in its indicated use, limitations, and possible complications of percutaneous needle techniques.
2. The introduction of the device into the body should be carried out under imaging control (ultrasound, X-Ray, CT, etc.). **NOTE: This product has not been tested for MR Imaging compatibility.**
3. Before using, inspect the device for damaged point, bent shaft or other imperfections that would prevent proper function. If the components are damaged or bent, DO NOT USE.

G. Potential Complications:

Potential complications are site specific and may consist of hematoma, hemorrhage, infection, adjacent tissue injury, pain; bleeding; hemoptysis; hemothorax; non-target tissue, organ or vessel perforation; and pneumothorax.

H. Equipment Required:

- Appropriate imaging modality
- Surgical gloves and drapes
- Local anesthetic
- Other equipment as necessary

I. Directions for Use:

1. Inspect the package and product for damage and expiration date. If undamaged and unexpired, open the package and transfer the product onto the sterile field utilizing aseptic technique.
2. Prior to use, determine the appropriate BARD® DUALOK® Localization Wire required for the specific localization to be performed. Using aseptic technique, remove the wire/needle cannula from the package. Before using, inspect the wire/needle cannula for damaged tips, bent shafts, deformed locking elements or other imperfections that would prevent proper function. If the wire/needle cannula are damaged, bent or do not move together easily, DO NOT USE.



Localization Procedure:

The localization procedure should be performed using appropriate aseptic techniques.

1. While holding the cannula hub stationary, grasp the proximal end of the wire and slowly withdraw the dual locking elements into the distal tip of the cannula. The wire is fully withdrawn when the band of tight twists on the wire are visible outside the cannula hub (See Figure 1).
2. Identify the desired tissue entry location and advance the cannula into the tissue. If needed, use the reference markings on the cannula to determine placement depth. Stop advancing when the distal tip of the cannula is just proximal to the desired location.
3. To set the wire locking elements, hold the needle cannula stationary while advancing the wire forward into the cannula hub. The locking elements are fully deployed when the band of tight twists on the wire are just inside the cannula locking hub (See Figure 2).

Repositioning or Removal of BARD® DUALOK® Localization Wire:

4. If the localization wire needs to be repositioned or removed, withdraw the locking elements into the needle cannula by holding the proximal end of the wire stationary and slowly advancing the needle cannula forward. The locking elements are fully withdrawn into the cannula when the bands of tight twists are visible outside the cannula hub. Then reposition the needle cannula to the desired location and re-set the locking element (See Step 3).
5. With the locking elements deployed, the localization wire is now locked into position. The cannula can now either be removed, or left in place to prevent accidental wire transection during surgery.

Figure 1**Figure 2**

Warranty

Bard Peripheral Vascular, Inc. warrants to the first purchaser of this product that this product will be free from defects in materials and workmanship for a period of one year from the date of first purchase and liability under this limited product warranty will be limited to repair or replacement of the defective product, in Bard Peripheral Vascular's sole discretion or refunding your net price paid. Wear and tear from normal use or defects resulting from misuse of this product are not covered by this limited warranty.

TO THE EXTENT ALLOWABLE BY APPLICABLE LAW, THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL BARD PERIPHERAL VASCULAR BE LIABLE TO YOU FOR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM YOUR HANDLING OR USE OF THIS PRODUCT.

Some states/countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your state/country.

An issue or revision date and revision number for these instructions are included for the user's information on the last page of this booklet. In the event 36 months have elapsed between this date and product use, the user should contact Bard Peripheral Vascular to see if additional product information is available.

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**Copyright ©2012 C. R. Bard Inc. All rights reserved.**
Copyright ©2012 C. R. Bard Inc. Tous droits réservés.
Copyright ©2012 C. R. Bard, Inc. Alle Rechte vorbehalten.
Copyright © 2012 C. R. Bard, Inc. Tutti i diritti riservati.
Copyright ©2012, C. R. Bard, Inc. Reservados todos los derechos.
Copyright © 2012 C.R. Bard, Inc. Alle rechten voorbehouden.