BARD® Blunt Tip Stiffening Cannula Instructions for Use

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BIOPSY SYSTEMS

Instructions for Use

Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

A. General Information and Device Description:

The BARD® Blunt Tip Stiffening Cannula consists of a cannula and hub.

B. How Supplied:

The BARD® Blunt Tip Stiffening Cannula is supplied sterile and non-pyrogenic unless the package has been damaged or opened. Sterilized using Ethylene Oxide. For single use only. Do not reuse. Do not resterilize.

Catalogue Number	Gauge Size and Needle Length
48500	20g x 5cm
48700	20g x 7cm
48900	20g x 9cm

C. Indications for Use:

This device is intended to be used as an aid for the GHIATAS® Beaded Breast Localization Wire or CHESBROUGH® Breast Localization Wire.

D. Contraindications:

This device is not intended for use except as indicated above.

E. Warnings:

- The BARD® Blunt Tip Stiffening Cannula has been designed for single use only. Reusing this medical device bears the risk of cross-patient contamination as medical devices – particularly those with long and small lumina, joints, and/or crevices between components – are difficult or impossible to clean once body fluids or tissues with potential pyrogenic or microbial contamination have had contact with the medical device for an indeterminable period of time. The residue of biological material can promote the contamination of the device with pyrogens or microorganisms which may lead to infectious complications.
- Do not resterilize the BARD® Blunt Tip Stiffening Cannula. After resterilization, the sterility of the product is not guaranteed because of an indeterminable degree of potential pyrogenic or microbial contamination which may lead to infectious complications. Cleaning, reprocessing and/or resterilization of the present medical device increases the probability that the device will malfunction due to potential adverse effects on components that are influenced by thermal and/or mechanical changes.

Note: After use, this product may be a potential biohazard. Handle and dispose of in accordance with acceptable medical practice and local, state, and federal laws and regulations.

F. Precautions:

- This device should only be used by a physician trained in its indicated use, limitations, and possible complications.
- The introduction of the device into the body should be carried out under imaging control (ultrasound or x-ray).
- Before using, inspect the device for damage that would prevent proper function. If the components are damaged or bent, DO NOT USE.

G. Potential Complications:

Potential complications associated with core biopsy procedures are site specific and include, but are not limited to: hematoma; hemorrhage; infection; adjacent tissue injury; pain; bleeding; hemoptysis; hemothorax; non-target tissue, organ or vessel perforation; and pneumothorax.

H. Directions for Use:

Prior to surgery the BARD® Blunt Tip Stiffening Cannula may be guided over the GHIATAS® Beaded Breast Localization Wire or CHESBROUGH® Breast Localization Wire to aid in providing a palpable indication of the course of the wire and location of the wire tip.

For more information refer to your GHIATAS® Beaded Breast Localization Wire or CHESBROUGH® Breast Localization Wire Instructions for Use.

Warranty:

Bard Peripheral Vascular warrants to the first purchaser of this product that this product will be free from defects in materials and workmanship for a period of one year from the date of first purchase and liability under this limited product warranty will be limited to repair or replacement of the defective product, in Bard Peripheral Vascular's sole discretion or refunding your net price paid. Wear and tear from normal use or defects resulting from misuse of this product are not covered by this limited warranty.

TO THE EXTENT ALLOWABLE BY APPLICABLE LAW, THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL BARD PERIPHERAL VASCULAR BE LIABLE TO YOU FOR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM YOUR HANDLING OR USE OF THIS PRODUCT.

Some states/countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your state/country.

An issue or revision date and a revision number for these instructions are included for the user's information on the last page of this booklet. In the event 36 months have elapsed between this date and product use, the user should contact Bard Peripheral Vascular to see if additional product information is available.

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Bard, Chesbrough, Ghiatas, and UltraWire are trademarks and/or registered trademarks of C. R. Bard, Inc. or an affiliate.
Copyright © 2012 C. R. Bard, Inc. All rights reserved.
Copyright © 2012 C. R. Bard, Inc. Tous droits réservés.
Copyright © 2012 C. R. Bard, Inc. Alle Rechte vorbehalten.
Copyright © 2012 C. R. Bard, Inc. Tutti i diritti riservati.
Copyright © 2012 C. R. Bard, Inc. Reservados todos los derechos.