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BIOPSY SYSTEMS


Instructions for Use

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

A. General Information and Device Description:

The CHESBROUGH® Breast Localization Wire consists of an introducer needle and a semi-rigid dual barb localization wire. The introducer needle is comprised of a plastic molded hub, 1cm depth reference marks, and an ultrasound enhancement on the distal end to aid in needle placement. The semi-rigid dual barb localization wire has a bead on the distal end to aid in visualization and a palpable deployment bead to signify the point at which the barbs at the distal end of the wire will be deployed.

Figure 1

1. Introducer needle
2. Ultrasound enhancement
3. 1cm depth reference marks
4. Orientation marks
5. Barbs
6. Bead
7. Palpable deployment bead
8. Localization wire

B. How Supplied:

The CHESBROUGH® Breast Localization Wire is supplied sterile and non-pyrogenic unless the package has been opened or damaged. Sterilized using Ethylene Oxide. **For single use only. Do not reuse. Do not resterilize.**

C. Indications for Use:

This device is intended for use during breast lesion surgery as a guide for the surgeon to follow in the excision of the lesion.

D. Contraindications for Use:

None known.

E. Warnings:

1. Once the barbs have been deployed into the breast, the wire must be removed surgically. Do not attempt to reposition, move, or pull on the wire or damage/breakage may result.
2. The wire is meant to be a guide only. DO NOT use it as a retractor.
3. Exercise caution during surgical excision of the lesion to avoid cutting the wire with a scalpel.
4. The wire should only be advanced through the needle from the hub to the needle point. Do not attempt to pull wire backwards into or inside the needle or barb damage may result.
5. Do not attempt to reshape the barbs of the wire in any manner; this may cause the barbs to fracture. If the barbs of the wire are deformed or bent incorrectly, discard product.
6. Caution should be exercised when using any localization wire on patients with breast prostheses so as not to puncture them during placement or transportation.
7. The CHESBROUGH® Breast Localization Wire has been designed for single use only. Reusing this medical device bears the risk of cross-patient contamination as medical devices – particularly those with long and small lumina, joints, and/or crevices between components – are difficult or impossible to clean once body fluids or tissues with potential pyrogenic or microbial contamination have had contact with the medical device for an indeterminable period of time. The residue of biological material can promote the contamination of the device with pyrogens or microorganisms which may lead to infectious complications.
8. Do not resterilize the CHESBROUGH® Breast Localization Wire. After resterilization, the sterility of the product is not guaranteed because of an indeterminable degree of potential pyrogenic or microbial contamination which may lead to infectious complications. Cleaning, reprocessing and/or resterilization of the present medical device increases the probability that the device will malfunction due to potential adverse effects on components that are influenced by thermal and/or mechanical changes.

NOTE: After use, this product may be a potential biohazard. Handle and dispose of in accordance with acceptable medical practice and applicable local, state, and federal laws and regulations.

F. Precautions:

1. This device should only be used by a physician trained in its indicated use, limitations, and possible complications of percutaneous needle techniques.
2. The introduction of the device into the body should be carried out under imaging control (ultrasound or x-ray).
3. Before using, inspect the device for damage that would prevent proper function. If the components are damaged or bent, DO NOT USE.

G. Potential Complications:

Potential complications are site specific and may consist of hematoma; hemorrhage; infection; adjacent tissue injury; pain; bleeding; hemoptysis; hemothorax; non-target tissue, organ or vessel perforation; and pneumothorax.

H. Equipment Required:


- Appropriate imaging modality
- Surgical gloves and drapes
- Local anesthetic
- Other equipment as necessary

I. Directions for Use:

1. Inspect the package and product for damage and expiration date. If undamaged and unexpired, open the package and transfer the product onto the sterile field utilizing aseptic technique.
2. Remove the wire from the introducer needle and verify that the product was not damaged during shipping (Figure 2).

NOTE: The localization wire is packaged in the needle with the barbs facing out of the needle hub to prevent damage to the barbs. At the time of placement, the localization wire will need to be turned around so that the barbs slide through the needle hub first.

Figure 2


3. Locate the lesion using the appropriate ultrasound or x-ray imaging technique.
4. Insert the localization needle into the breast, directing it to the lesion (preferably parallel to the chest wall to reduce the possible risk of pneumothorax). Use the depth reference marks to position the needle point into the breast lesion (depth reference marks are 1cm apart) (Figure 3).

Figure 3


5. Confirm needle placement with x-ray and/or ultrasound.
6. If necessary, reposition the needle and reconfirm placement.
7. The localization wire may then be inserted into the introducer needle by sliding the wire into the needle hub being careful not to snag the proximal barb outside the hub. Advance the wire until the beginning of the palpable deployment bead on the wire is at the proximal end of the needle hub. To insure the barbs will deploy in the proper direction, line the barbs up with the orientation marks on the hub when feeding the wire into the needle (Figure 4).

Figure 4

8. To deploy the barbs, the localization wire should be held in place and the needle withdrawn approximately 1.5cm (Figure 5). **NOTE:** The barbs will deploy when the palpable deployment bead on the wire is completely inside the needle hub. Any effort to reposition the wire after the barbs are deployed could damage/break the wire.

Figure 5

9. The needle can now be removed and an x-ray and/or ultrasound can be taken to confirm the placement of the barbs. **NOTE:** The bead on the wire offers a visual and palpable reference to the location of the barbs and lesion.
10. The remaining exposed wire should be secured to the skin surface, using an appropriate method to prevent the wire from moving during transport to surgery (Figure 6). The wire clip may be placed on the wire at the skin surface to help avoid wire migration. Use sufficient pressure to place the wire clip on the wire (Figure 7).

Figure 6



Figure 7

11. Prior to surgery the BARD® Blunt Tip Stiffening Cannula may be guided over the CHESBROUGH® Breast Localization Wire to aid in providing a palpable indication of the course of the wire and location of the wire tip.

Warranty

Bard Peripheral Vascular, Inc. warrants to the first purchaser of this product that this product will be free from defects in materials and workmanship for a period of one year from the date of first purchase and liability under this limited product warranty will be limited to repair or replacement of the defective product, in Bard Peripheral Vascular's sole discretion or refunding your net price paid. Wear and tear from normal use or defects resulting from misuse of this product are not covered by this limited warranty.

TO THE EXTENT ALLOWABLE BY APPLICABLE LAW, THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL BARD PERIPHERAL VASCULAR BE LIABLE TO YOU FOR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM YOUR HANDLING OR USE OF THIS PRODUCT.

Some states/countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your state/country.

An issue or revision date and revision number for these instructions are included for the user's information on the last page of this booklet. In the event 36 months have elapsed between this date and product use, the user should contact Bard Peripheral Vascular to see if additional product information is available.

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Bard, Chesbrough, Ghiatas ve UltraWire are trademarks and/or registered trademarks of C. R. Bard, Inc. or its affiliate. All other trademarks are the property of their respective owners.

Bard, Chesbrough, Ghiatas and UltraWire are trademarks and/or registered trademarks of C. R. Bard, Inc. or its affiliate. All other trademarks are the property of their respective owners.

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Copyright © 2012 C. R. Bard, Inc. All rights reserved.

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Bard, Chesbrough, Ghiatas, and UltraWire are trademarks and/or registered trademarks of C. R. Bard, Inc. or an affiliate. All other trademarks are the property of their respective owners.
Copyright © 2012 C. R. Bard, Inc. Tous droits réservés.
Copyright © 2012 C. R. Bard, Inc. Alle Rechte vorbehalten.
Copyright © 2012 C. R. Bard, Inc. Tutti i diritti riservati.
Copyright © 2012, C. R. Bard, Inc. Reservados todos los derechos.
Copyright © 2012 C. R. Bard, Inc. Alle rechten voorbehouden.