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
BIOPSY SYSTEMS

Instructions for Use

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

Figure 1

1. Breast localization device
2. Ultrasound enhancement
3. 1cm depth reference marks
4. Trigger
5. Barbs
6. Bead
7. Localization wire

A. General Information and Device Description:

The ULTRA WIRE[®] One-Handed Breast Localization Device consists of a plastic molded handle, an introducer needle with 1cm depth reference marks, and an ultrasound enhancement on the distal end to aid in needle placement. A semi-rigid dual barb localization wire is loaded inside the handle. The handle has a trigger which when activated will cause the needle to retract back into the handle exposing the dual barb localization wire. The wire contains a visual bead between the barbs.

B. How Supplied:

The ULTRA WIRE[®] One-Handed Breast Localization Device is supplied sterile and non-pyrogenic unless the package has been damaged or opened. Sterilized using Ethylene Oxide. **For single use only. Do not reuse. Do not resterilize.**

C. Indications for Use:

This device is intended for use during breast lesion surgery as a guide for the surgeon to follow in the excision of the lesion.

D. Contraindications for Use:

None known.

E. Warnings:

1. **Once the barbs have been deployed into the breast, the wire must be removed surgically. Do not attempt to reposition, move, or pull on the wire or damage/breakage may result.**
2. **The wire is meant to be a guide only. DO NOT use it as a retractor.**
3. **Exercise caution during surgical excision of the lesion to avoid cutting the wire with a scalpel.**
4. **Caution should be exercised when using any localization wire on patients with breast prostheses so as not to puncture them during placement or transportation.**
5. **The ULTRA WIRE[®] One-Handed Breast Localization Device has been designed for single use only. Reusing this medical device bears the risk of cross-patient contamination as medical devices – particularly those with long and small lumina, joints, and/or crevices between components – are difficult or impossible to clean once body fluids or tissues with potential pyrogenic or microbial contamination have had contact with the medical device for an indeterminable period of time. The residue of biological material can promote the contamination of the device with pyrogens or microorganisms which may lead to infectious complications.**
6. **Do not resterilize the ULTRA WIRE[®] One-Handed Breast Localization Device. After resterilization, the sterility of the product is not guaranteed because of an indeterminable degree of potential pyrogenic or microbial contamination which may lead to infectious complications. Cleaning, reprocessing and/or resterilization of the present medical device increases the probability that the device will malfunction due to potential adverse effects on components that are influenced by thermal and/or mechanical changes.**

NOTE: After use, this product may be a potential biohazard. Handle and dispose of in accordance with acceptable medical practice and local, state, and federal laws and regulations.

F. Precautions:

1. This device should only be used by a physician trained in its indicated use, limitations, and possible complications of percutaneous needle techniques.
2. The introduction of the device into the body should be carried out under imaging control (ultrasound or x-ray).
3. Before using, inspect the device for damage that would prevent proper function. If the components are damaged or bent, **DO NOT USE.**

G. Potential Complications:

Potential complications are site specific and may consist of hematoma, hemorrhage, infection, adjacent tissue injury, pain, bleeding, hemoptysis, hemothorax, non-target tissue, organ or vessel perforation, and pneumothorax.

H. Equipment Required:

- Appropriate imaging modality
- Surgical gloves and drapes
- Local anesthetic
- Other equipment as necessary

I. Directions for Use:

1. Inspect the package and product for damage and expiration date. If undamaged and unexpired, open the package and transfer the product onto the sterile field utilizing aseptic technique.
2. Locate the lesion using the appropriate ultrasound or x-ray imaging technique.
3. Remove tip protector from needle by pulling straight off to prevent tip damage.
4. Insert the localization needle into the breast, directing it to the lesion (preferably parallel to the chest wall to reduce the possible risk of pneumothorax). Use the depth reference marks to position the needle point into the breast lesion (depth reference marks are 1cm apart).
5. Confirm needle placement with ultrasound and/or x-ray.
6. If necessary, reposition the needle and reconfirm placement.
7. Deploy the wire by pulling the trigger slightly rearward (away from the needle) and depress the trigger (**Figure 2**). The needle will retract into the handle automatically (**Figure 3**). **NOTE:** The needle must be straight in order to retract. Bending the needle may prevent the device from functioning properly.

Figure 2



Figure 3

8. Confirm placement of the wire with imaging.
9. The remaining exposed wire should be secured to the skin surface, using an appropriate method to prevent the wire from moving during transport to surgery.

Warranty

Bard Peripheral Vascular, Inc. warrants to the first purchaser of this product that this product will be free from defects in materials and workmanship for a period of one year from the date of first purchase and liability under this limited product warranty will be limited to repair or replacement of the defective product, in Bard Peripheral Vascular's sole discretion or refunding your net price paid. Wear and tear from normal use or defects resulting from misuse of this product are not covered by this limited warranty.

TO THE EXTENT ALLOWABLE BY APPLICABLE LAW, THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL BARD PERIPHERAL VASCULAR BE LIABLE TO YOU FOR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM YOUR HANDLING OR USE OF THIS PRODUCT.


Some states/countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your state/country.

An issue or revision date and a revision number for these instructions are included for the user's information on the last page of this booklet. In the event 36 months have elapsed between this date and product use, the user should contact Bard Peripheral Vascular to see if additional product information is available.

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Bard, Chesbrough, Ghiatas, and UltraWire are trademarks and/or registered trademarks of C. R. Bard, Inc. or an affiliate.
Copyright © 2012 C. R. Bard, Inc. All rights reserved.
Copyright © 2012 C. R. Bard, Inc. Tous droits réservés.
Copyright © 2012 C. R. Bard, Inc. Alle Rechte vorbehalten.
Copyright © 2012 C. R. Bard, Inc. Tutti i diritti riservati.
Copyright © 2012 C. R. Bard, Inc. Todos los derechos reservados.