BARD® MAX-CORE®
Disposable Core Biopsy Instrument
Instructions for Use

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BIOPSY SYSTEMS

Instructions for Use:

Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

A. General Information and Device Description:

The BARD® MAX-CORE® Biopsy Instrument is a single use core biopsy device. It is available in several needle gauge sizes and lengths. The side and rear actuator buttons are color coded according to the various gauge sizes, e.g., Yellow=20 gauge, Pink=18 gauge, Purple=16 gauge and Green=14 gauge.

Catalogue Number	Gauge Size and Needle Length	Length of Sample Notch	Penetration Depth
MC1410	14g (2.1mm) x 10cm (100mm)	1.9cm (19mm)	22mm
MC1416	14g (2.1mm) x 16cm (160mm)	1.9cm (19mm)	22mm
MC1610	16g (1.7mm) x 10cm (100mm)	1.9cm (19mm)	22mm
MC1616	16g (1.7mm) x 16cm (160mm)	1.9cm (19mm)	22mm
MC1810	18g (1.2mm) x 10cm (100mm)	1.8cm (18mm)	22mm
MC1816	18g (1.2mm) x 16cm (160mm)	1.8cm (18mm)	22mm
MC1820	18g (1.2mm) x 20cm (200mm)	1.8cm (18mm)	22mm
MC1825	18g (1.2mm) x 25cm (250mm)	1.8cm (18mm)	22mm
MC2010	20g (0.9mm) x 10cm (100mm)	1.8cm (18mm)	22mm
MC2016	20g (0.9mm) x 16cm (160mm)	1.8cm (18mm)	22mm
MC2020	20g (0.9mm) x 20cm (200mm)	1.8cm (18mm)	22mm

B. How Supplied:

The product is supplied sterile and non-pyrogenic unless the package has been opened or damaged. **Sterilized using Ethylene Oxide. For single use only. Do Not Reuse. Do Not Resterilize.**

C. Indications for Use:

The core needle biopsy device is intended for use in obtaining biopsies from soft tissues such as liver, kidney, prostate, spleen, lymph nodes and various soft tissue tumors. It is not intended for use in bone.

D. Contraindications:

Good medical judgment should be exercised in considering biopsy on patients who are receiving anticoagulant therapy or who have a bleeding problem.

E. Warnings:

1. Post-biopsy patient care may vary with the biopsy technique utilized and the individual patient's physiological condition. Observation of vital signs and other precautions should be taken to avoid and/or treat potential complications that may be associated with biopsy procedures.
2. The collection of multiple needle cores may help to ensure the detection of any cancer tissue. A "negative" biopsy in the presence of suspicious radiographic findings does not preclude the presence of carcinoma.
3. The BARD® MAX-CORE® Biopsy Instrument has been designed for single use only. Reusing this medical device bears the risk of cross-patient contamination as medical devices – particularly those with long and small lumina, joints, and/or crevices between components – are difficult or impossible to clean once body fluids or tissues with potential pyrogenic or microbial contamination have had contact with the medical device for an indeterminable period of time. The residue of biological material can promote the contamination of the device with pyrogens or microorganisms which may lead to infectious complications.
4. Do not resterilize the BARD® MAX-CORE® Biopsy Instrument. After resterilization, the sterility of the product is not guaranteed because of an indeterminable degree of potential pyrogenic or microbial contamination which may lead to infectious complications. Cleaning, reprocessing and/or resterilization of the present medical device increases the probability that the device will malfunction due to potential adverse effects on components that are influenced by thermal and/or mechanical changes.

Note: If collecting multiple samples, inspect the needle for damaged point, bent shaft or other imperfections after each sample is collected. Do not use needle if any imperfection is noted.

Note: After use, this product may be a potential biohazard. Handle and dispose of in accordance with acceptable medical practice and with applicable local, state, and federal laws and regulations.

F. Precautions:

1. This product should be used by a physician who is completely familiar with the indications, contraindications, limitations, typical findings and possible side effects of core needle biopsy, in particular, those relating to the specific organ being biopsied.
2. The introduction of the needle into the body should be carried out under imaging control (ultrasound, X-Ray, CT, etc.).
3. Never test the product by firing into the air. Damage may occur to the needle/cannula tip and could result in patient and/or user injury.
4. Unusual force applied to the stylet or unusual resistance against the stylet while extended out of the supportive cannula may cause the stylet to bend at the specimen notch. A bent specimen notch may interfere with the needle function.

G. Potential Complications:

Potential complications associated with core biopsy procedures are site specific and include, but are not limited to: hematoma; hemorrhage; infection; adjacent tissue injury; pain; bleeding; hemoptysis; hemothorax; non-target tissue, organ or vessel perforation; and air embolism. Air embolism is a rare but serious potential complication of lung biopsy procedures. Rapid deterioration of neurological status and/or cardiac arrhythmia may be indicative of air embolism. Prompt diagnosis and treatment must be considered if the patient exhibits signs or symptoms of air embolism.

H. Equipment Required:

- Appropriate imaging modality accessories
- Surgical gloves and drapes
- Local anesthetic
- BARD® TRUGUIDE® Coaxial Cannula (optional)
- Scalpel
- Sample collection container
- Other equipment as necessary

I. Directions for Use:**BARD® MAX-CORE® Biopsy Instrument preparation:**

Before using, inspect the needle for damaged point, bent shaft or other imperfections that would prevent proper function. If the needle is damaged or bent, DO NOT USE.

1. Using aseptic technique, remove the Instrument from its package. See Figure 1.

Note: To remove the protective needle sheath and yellow guard, you must energize (cock) Instrument.

2. Energize (cock) Instrument by pulling back on the top slide to withdraw the cannula and lock in place. See Figure 2. Then pull back on the bottom slide to withdraw the stylet and lock in place. See Figure 3. Remove protective needle sheath and yellow guard. Instrument is ready to fire when both slides are locked back.

Recommendation: For ease of insertion, puncture the skin with a scalpel at the entry site.

Biopsy Procedure:

The biopsy procedure must be performed using appropriate aseptic techniques.

1. Prepare site as required. Adequate anesthesia should be administered prior to incision of the skin.
2. Verify Instrument is energized (cocked). See Figure 3.

Note: Do not place fingers in front of cocking slides once Instrument is energized (cocked). Impeding cocking slides' movement will impact functionality.


3. Insert tip of needle to the point to be biopsied.
4. While maintaining Instrument's position and the needle orientation, depress the rear actuator button, or push the side actuator forward (direction of arrow), to cause both stylet and cannula to automatically advance.
5. Remove needle from patient and pull back on the top slide to withdraw the cannula and expose the biopsy specimen (See Figure 2). Remove the specimen.
6. If additional biopsies are required, pull back on the bottom slide to withdraw the stylet and repeat the procedure.

Figure 1

- As Packaged (Protective Needle Sheath and Yellow Guard Not Shown)
- Instrument **Not** Energized (Cocked)
- **Not** Ready to Fire

Figure 2

- Top Slide Locked Back
- Biopsy Sample Notch Exposed

Figure 3

- Top Slide and Bottom Slide Locked Back
- Instrument Energized (Cocked)
- Ready to Fire

Warranty

Bard Peripheral Vascular warrants to the first purchaser of this product that this product will be free from defects in materials and workmanship for a period of one year from the date of first purchase and liability under this limited product warranty will be limited to repair or replacement of the defective product, in Bard Peripheral Vascular's sole discretion or refunding your net price paid. Wear and tear from normal use or defects resulting from misuse of this product are not covered by this limited warranty.

TO THE EXTENT ALLOWABLE BY APPLICABLE LAW, THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL BARD PERIPHERAL VASCULAR BE LIABLE TO YOU FOR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM YOUR HANDLING OR USE OF THIS PRODUCT.

Some states/countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your state/country.

An issue or revision date and revision number for these instructions are included for the user's information on the last page of this booklet. In the event 36 months have elapsed between this date and product use, the user should contact Bard Peripheral Vascular to see if additional product information is available.

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**Bard, Max-Core, and TruGuide are trademarks and/or registered trademarks of C. R. Bard, Inc. or an affiliate.**
Bard, Max-Core and TruGuide are C. R. Bard, Inc. or its affiliated organization's trademarks and/or registered trademarks.
Copyright ©2012 C. R. Bard Inc. All rights reserved.
Copyright ©2012 C. R. Bard Inc. Tous droits réservés.
Copyright ©2012 C. R. Bard Inc. Alle Rechte vorbehalten.
Copyright ©2012 C. R. Bard Inc. Tutti i diritti riservati.