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
BARD® Needle Guide for use with the Lorad™ StereoGuide™ System

Instructions for Use

Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

A. General Information and Device Description:

The BARD® Needle Guide is used with the Lorad™ prone stereotactic biopsy table. The needle guides are a disposable, single patient use product.

B. How Supplied:

The BARD® Needle Guide for use with the Lorad™ StereoGuide™ System is supplied sterile unless the package has been damaged or opened. Sterilized using Ethylene Oxide. **For single use only. Do not reuse. Do not resterilize.**

C. Indications for Use:

The BARD® Needle Guide is intended for use as a needle guide during stereotactic guided breast biopsy procedures.

D. Contraindications:

None known.

E. Warnings:

1. The BARD® Needle Guide has been designed for single use only. Reusing this medical device bears the risk of cross-patient contamination as medical devices – particularly those with long and small lumina, joints, and/or crevices between components – are difficult or impossible to clean once body fluids or tissues with potential pyrogenic or microbial contamination have had contact with the medical device for an indeterminable period of time. The residue of biological material can promote the contamination of the device with pyrogens or microorganisms which may lead to infectious complications.
2. Do not resterilize the BARD® Needle Guide. After resterilization, the sterility of the product is not guaranteed because of an indeterminable degree of potential pyrogenic or microbial contamination which may lead to infectious complications. Cleaning, reprocessing and/or resterilization of the present medical device increases the probability that the device will malfunction due to potential adverse effects on components that are influenced by thermal and/or mechanical changes.

Note: After use, this product may be a potential biohazard. Handle and dispose of in accordance with acceptable medical practice and applicable local, state, and federal laws and regulations.

Warranty

Bard Peripheral Vascular warrants to the first purchaser of this product that this product will be free from defects in materials and workmanship for a period of one year from the date of first purchase and liability under this limited product warranty will be limited to repair or replacement of the defective product, in Bard Peripheral Vascular's sole discretion or refunding your net price paid. Wear and tear from normal use or defects resulting from misuse of this product are not covered by this limited warranty.

TO THE EXTENT ALLOWABLE BY APPLICABLE LAW, THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL BARD PERIPHERAL VASCULAR BE LIABLE TO YOU FOR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM YOUR HANDLING OR USE OF THIS PRODUCT.

Some states/countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your state/country.

An issue or revision date and a revision number for these instructions are included for the user's information on the last page of this booklet. In the event 36 months have elapsed between this date and product use, the user should contact Bard Peripheral Vascular, Inc. to see if additional product information is available.

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Bard is a trademark and/or registered trademark of C. R. Bard, Inc. or an affiliate. All other trademarks are the property of their respective owners.
Copyright © 2012 C. R. Bard, Inc. All rights reserved.
Copyright © 2012 C. R. Bard, Inc. Tous droits réservés.
Copyright © 2012 C. R. Bard, Inc. Alle Rechte vorbehalten.
Copyright © 2012 C. R. Bard, Inc. Tutti i diritti riservati.
Copyright © 2012 C. R. Bard, Inc. Todos los derechos reservados.