ENGLISH

REF GMU11T
 MK2011

Instructions For Use

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.


Figure 1 – Biopsy Site Marker

Description

The GEL MARK ULTRA™ Biopsy Site Marker is a sterile, disposable applicator containing resorbable polylactic acid/polyglycolic acid (PLA/PGA) pellets in the following configuration:

- 11 resorbable pellets,
- One with embedded radiopaque wireform in the sixth position.
 MK2011: 316L stainless steel wireform with an omega shape
 GMU11T: Titanium wireform in an “S” shape

Product Code	Wireform Material	Wireform Shape	Probe Compatibility
GMU11T	Titanium		11G Mammotome™
MK2011	316L stainless steel		11G Mammotome™

The wireform is intended for long-term radiographic marking of the biopsy site. The pellets are visible via ultrasound for approximately 6 weeks and are essentially resorbed in approximately 12 weeks.

The syringe-like applicator fits within the Mammotome™ to access the biopsy cavity. Pellets are deployed from the applicator through the biopsy device into the biopsy cavity.

Intended Use

The GEL MARK ULTRA™ is intended to radiographically mark breast tissue during a percutaneous breast biopsy procedure.

Contraindications

None known

Warnings

- **This device has been designed for single use only. Reusing this medical device bears the risk of cross-patient contamination as medical devices – particularly those with long and small lumina, joints, and/or crevices between components – are difficult or impossible to clean once body fluids or tissues with potential pyrogenic or microbial contamination have had contact with the medical device for an indeterminable period of time. The residue of biological material can promote the contamination of the device with pyrogens or microorganisms which may lead to infectious complications.**
- **Do not resterilize. After resterilization, the sterility of the product is not guaranteed because of an indeterminable degree of potential pyrogenic or microbial contamination which may lead to infectious complications. Cleaning, reprocessing and/or resterilization of the present medical device increases the probability that the device will malfunction due to potential adverse effects on components that are influenced by thermal and/or mechanical changes.**
- **Patients with a known hypersensitivity to the materials listed in the device description may suffer an allergic reaction to this implant.**
- **The GEL MARK ULTRA™ is not recommended for use in patients with breast implants.**
- **Do not use in the presence of infection.**
- **Avoid the use of excessive force during removal of the applicator to prevent breakage of the applicator tip.**

Precautions

- The device should only be used by physicians trained in percutaneous biopsy procedures.
- Do not use this product if the sterile barrier has been previously opened or if the package is damaged.

- This product is provided sterile, and is intended for single use only. Do Not Resterilize.
- Do not kink the flexible tube.
- Maintain correct alignment of the yellow indicator key (C in Figure 1) with the white arrow of the biopsy probe when dispensing pellets.
- Ensure that all pellets are dispensed.
- For model MK2011, with stainless steel marker:
 Published studies with comparably sized 316L stainless steel biopsy markers have shown no marker movement and insignificant marker heating when tested in an MRI System with a 1.5T field strength.

COMPLICATIONS

Potential complications that may be associated with the use of the GEL MARK ULTRA™ are similar to those associated with the use of other biopsy marking devices.

DIRECTIONS FOR USE

1. Perform a “dry tap” to ensure that the sample notch is clear of tissue.
2. Using standard aseptic technique, remove the applicator from the package and inspect for damage. Remove tip cover.
3. Insert the applicator into the biopsy device. Advance the applicator to align the yellow indicator key (C) with the white arrow on the thumbwheel of the biopsy probe. Refer to Figure 2.


Figure 2 – Applicator alignment in Probe

4. Fully seat the applicator in the biopsy probe. Verify proper alignment of the applicator aperture with the biopsy chamber by ensuring proper alignment of the yellow indicator key (C) with the white arrow on the thumbwheel of the biopsy probe (Figure 2).
 5. Immediately deploy all pellets by advancing the plunger (A).
- Note:** Avoid pulling the applicator handle during deployment which could cause deployment failure.
6. Rotate the biopsy device cannula 180° to position the sample chamber away from the dispensed pellets.
 7. Remove the applicator from the biopsy device.
- Caution:** If resistance is felt while removing the applicator, remove the entire probe/applicator assembly from the patient. Failure to do so may result in breakage of the applicator tip.
8. Close the sample chamber and remove the biopsy device, per manufacturer’s instructions.
 9. Dispose of the applicator properly.
 10. Confirm final marker position with imaging.

WARRANTY

Bard Peripheral Vascular warrants to the first purchaser of this product that this product will be free from defects in materials and workmanship for a period of one year from the date of first purchase and liability under this limited product warranty will be limited to repair or replacement of the defective product, in Bard Peripheral Vascular’s sole discretion or refunding your net price paid. Wear and tear from normal use or defects resulting from misuse of this product are not covered by this limited warranty.

TO THE EXTENT ALLOWABLE BY APPLICABLE LAW, THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL BARD PERIPHERAL VASCULAR BE LIABLE TO YOU FOR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM YOUR HANDLING OR USE OF THIS PRODUCT.

Some states/countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your state/country.

An issue or revision date and a revision number for these instructions are included for the user’s information on the last page of this booklet. In the event 36 months have elapsed between this date and product use, the user should contact Bard Peripheral Vascular, Inc. to see if additional product information is available.

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

**Do Not Use if Temperature Indicator is Black**

Ne pas utiliser si l'indicateur thermique est noir.
 Bei schwarz gefärbter Temperaturanzeige nicht verwenden
 Non utilizzare se l'indicatore della temperatura è nero.
 No usar si el punto de temperatura tiene color negro.
 Niet gebruiken wanneer de temperatuursindicator zwart is.
 Não utilize o indicador de temperatura estiver preto
 Μη χρησιμοποιείτε εάν ο δείκτης θερμοκρασίας είναι μαύρος.
 Må ikke anvendes, hvis temperaturpletten er sort
 Använd inte produkten om temperaturindikatorn är svart.
 Älä käyttää, jos lämpötilan ilmaisin on musta
 Skäl ikke brukes hvis temperaturindikatoren er svart.
 Nie używać, jeśli wskaźnik temperatury jest czarny
 Ne használnia, ha a hőmérséklet-kijelző fekete!
 Produkt nepoužívejte, pokud je indikátor teploty černý
 Sıcaklık Göstergesi Siyah Renkliyse Kullanmayın
 Запрещается использовать, если индикатор температуры имеет черный цвет.
 如果溫度指示標籤是黑色的，請勿使用

†

**Gauge Size**

Diámetro
 Gauge-Maß
 Dimensioni gauge
 Tamanho de calibre
 Naalddikte
 Calibre
 Μέγεθος διαμετρήματος
 Gauge-størrelse
 Gauge-storlek
 G-koko
 Målerstørrelse
 Rozmiar igły
 Gauge méret
 Velikost měřky
 Kalibre Boyutu
 Калибр
 直徑尺寸

**Upper Limit of Temperature**

Limite supérieure de température
 Obere Temperaturgrenze
 Temperatura massima
 Límite superior de temperatura
 Bovengrens temperatuur
 Limite superior de temperatura
 Ανώτατο όριο θερμοκρασίας
 Øvre temperaturgrense
 Högsta temperatur
 Lämpötilan yläraja
 Øvre temperaturgrense
 Maksymalna dopuszczalna temperatura
 A hőmérséklet felső határértéke
 Horní mez teploty
 Üst Sıcaklık Limiti
 Верхний предел температуры
 溫度上限

**Not Made with Natural Rubber Latex**

Fabriqué sans latex naturel
 Nicht aus Naturgummilatem hergestellt.
 Non prodotto con lattice di gomma naturale
 Este producto no se fabrica con látex de caucho natural
 Niet gemaakt met latex van natuurlijk rubber
 Não contém látex de borracha natural
 Δεν κατασκευάζεται από φυσικό ελαστικό λάτεξ
 Ikke fremstillet med naturlig gummilatem
 Ej tillverkad med naturgummilatem
 Valmistuksessa ei ole käytetty luonnonkumilateksia
 Dette produktet er ikke produsert med naturlig gummilatem
 Produkt nie jest wytwarzany z lateksu kauczuku naturalnego
 Nem tartalmaz természetes gumilateket.
 Nemí vyroben z prírodného pryžového latexu
 Doğal Kauçuk Lateksten Üretilmemiştir
 Не содержит натуральный каучуковый латекс
 製造未採用天然膠乳

**Sterilized Using Irradiation**

Stérilisé par irradiation
 Durch Bestrahlung sterilisiert
 Sterilizzato mediante irradiazione
 Esterilizado mediante irradiación
 Gesteriliseerd met behulp van straling
 Esterilizado por Irradiação
 Αποστειρωμένο με χρήση ακτινοβολίας
 Steriliseret ved brug af bestråling
 Steriliserad med strålning
 Steriloitu säteilyttämällä
 Steriliseret ved stråling
 Produkt sterylizowany przez napromienianie
 Sugárzással sterilizált
 Sterilizováno ozářením
 Radyoaktif Işınla Sterilize Edilmiştir
 Стерилизовано радиацией
 已使用放射線消毒

**Do Not Use If Package Is Damaged**

Ne pas utiliser si l'emballage est endommagé
 Nicht verwenden, wenn die Verpackung beschädigt ist
 Non utilizzare se la confezione è danneggiata
 No usar si el envase está dañado
 Niet gebruiken wanneer de verpakking beschadigd is
 Não utilizar se a embalagem estiver danificada
 Μη χρησιμοποιείτε το προϊόν εάν η συσκευασία του έχει υποστεί ζημιά
 Må ikke bruges, hvis emballagen er beskadiget
 Använd inte produkten om förpackningen skadats
 Ei saa käyttää, jos pakkaus on vaurioitunut
 Skäl ikke brukes hvis pakningen er skadet
 Nie stosować, jeśli opakowanie jest uszkodzone
 Ne használnia fel, ha a csomagolás sérült
 Nepoužívejte, pokud je obal poškozen
 Paket Hasarlıysa Kullanmayınız
 Не использовать, если упаковка повреждена
 如果包裝已受損，請勿使用

**Authorised Representative in the European Community**

Représentant autorisé au sein de la Communauté européenne
 Berechtigter Vertreter in der Europäischen Gemeinschaft
 Rappresentante autorizzato nella Comunità Europea
 Representante autorizado en la Comunidad Europea
 Gemachtigde vertegenwoordiger binnen de Europese Gemeenschap
 Representante autorizado na Comunidade Europeia
 Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα
 Autoriseret repræsentant i Det Europæiske Fællesskab
 Auktoriserad representant i EG
 Valtuutettu edustaja Euroopan yhteisössä
 Autoriseret representant i EU
 Autoryzowany przedstawiciel na terenie Unii Europejskiej
 Bejegyzett képviselő az Európai Közösségben
 Autorizovaný zástupce pro Evropskou unii
 Avrurpa Topluluğu Yetkili Temsilcisi
 Уполномоченный представитель в Европейском сообществе
 歐洲共同體的授權代表

**Breast Tissue Marker**

Marqueur de tissu mammaire
 Brustgewebemarker
 Marcatore tessutale mammario
 Marcador de tejido mamario
 Mammaweefselmarkering
 Marcador de tecido mamário
 Δείκτης μαστικού ιστού
 Brystvævsmarkør
 Bröstvävnadsmarkör
 Rintakudosmarkkeri
 Brystvevmarkør
 Marker tkanki piersiowej
 Emlőszövet marker
 Značkováč prsní tkáň
 Meme Dokusu İşaretleyicisi
 Биологический маркер молочной железы
 乳房組織標記

Copyright © 2015 C.R. Bard, Inc. All rights reserved.

Copyright © 2015 C.R. Bard, Inc. Tous droits réservés.
Copyright © 2015 C.R. Bard, Inc. Alle Rechte vorbehalten.
Copyright © 2015 C.R. Bard, Inc. Tutti i diritti riservati.
Copyright © 2015 C.R. Bard, Inc. Todos los derechos reservados.
Copyright © 2015 C.R. Bard, Inc. Alle rechten voorbehouden.
Copyright © 2015 C.R. Bard, Inc. Todos os direitos reservados.
Πνευματικά δικαιώματα © 2015 C.R. Bard, Inc. Με την επιφύλαξη παντός δικαιώματος.
Copyright © 2015 C.R. Bard, Inc. Alle rettigheder forbeholdt.
Copyright © 2015 C.R. Bard, Inc. Eftertryck förbjudes.
Copyright © 2015 C.R. Bard, Inc. Kaikki oikeudet pidätetään.
Copyright © 2015 C.R. Bard, Inc. Med enerrett.
Copyright © 2015 C.R. Bard, Inc. Wszelkie prawa zastrzeżone.
Copyright © 2015 C.R. Bard, Inc. Minden jog fenntartva.
Copyright © 2015 C.R. Bard, Inc. Všechna práva vyhrazena.
Telif Hakki © 2015 C.R. Bard, Inc. Tüm hakları saklıdır.
Авторские права © 2015 Компания С. R. Bard, Inc. Все права защищены.
Copyright © 2015 C.R. Bard, Inc. 保留所有權利。
© 2015 C.R. Bard, Inc.

Bard and Gel Mark Ultra are trademarks and/or registered trademarks of C. R. Bard, Inc. or an affiliate. All other trademarks are the property of their respective owners.

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**Authorised Representative in
the European Community**
Bard Limited
Forest House
Tilgate Forest Business Park
Brighton Road, Crawley
West Sussex
RH11 9BP, UK

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Bard Gel Mark Ultra C. R. Bard, Inc.
SenoRx, Inc.