SENO[®]MARK[®] ULTRA

SENO[®]MARK[®] ULTRA Breast Tissue Marker for Use with ENCOR[®] Biopsy Probe

INSTRUCTIONS FOR USE

CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

DESCRIPTION

The SENOMARK[®] ULTRA Breast Tissue Marker consists of a disposable applicator containing:

- 3 resorbable polyglycolic acid (PGA) pads that are essentially resorbed in approximately 12 weeks. The center pad contains a Titanium or BioDur™ 108 wireform with an interwoven polyvinyl alcohol (PVA) polymer for permanent ultrasound enhancement. Please note the PVA polymer is not resorbed.
- 1 polyethylene glycol (PEG) push pellet.

The wireform is intended for long-term radiographic marking of the biopsy site. The PVA polymer is intended for long-term sonographic marking of the biopsy site.

The syringe-like applicator fits within the ENCOR[®] Biopsy Probe to access the biopsy cavity. Pads are deployed from the applicator through the biopsy probe and into the biopsy cavity.

| Product Code | Wireform Material | Wireform Shape | Probe Compatibility |
|--------------|----------------------|---|---------------------|
| SMEC7R | Titanium with PVA | Ribbon  | 7G |
| SMEC7C | BioDur™ 108 with PVA | Coil  | 7G |
| SMEC10R | Titanium with PVA | Ribbon  | 10G |
| SMEC10C | BioDur™ 108 with PVA | Coil  | 10G |
| SMEC12R | Titanium with PVA | Ribbon  | 12G |
| SMEC12C | BioDur™ 108 with PVA | Coil  | 12G |

INDICATIONS FOR USE

The SENOMARK[®] ULTRA Breast Tissue Marker is intended to radiographically and sonographically mark breast tissue during a percutaneous breast biopsy procedure.

CONTRAINDICATIONS

Patients with a known hypersensitivity to the materials listed in the device description may suffer an allergic reaction to this implant.

WARNINGS

1. Use caution when placing near a breast implant to avoid puncture of the implant capsule.
2. Avoid the use of excessive force during removal of the applicator to prevent breakage of the applicator tip.
3. Do not use in the presence of infection.
4. This device has been designed for single use only. Reusing this medical device bears the risk of cross-patient contamination as medical devices – particularly those with long and small lumina, joints, and/or crevices between components – are difficult or impossible to clean once body fluids or tissues with potential pyrogenic or microbial contamination have had contact with the medical device for an indeterminable period of time. The residue of biological material can promote the contamination of the device with pyrogens or microorganisms which may lead to infectious complications.

5. Do not resterilize. After resterilization, the sterility of the product is not guaranteed because of an indeterminable degree of potential pyrogenic or microbial contamination which may lead to infectious complications. Cleaning, reprocessing and/or resterilization of the present medical device increases the probability that the device will malfunction due to potential adverse effects on components that are influenced by thermal and/or mechanical changes.

PRECAUTIONS

1. This product should only be used by a physician who is completely familiar with the indications, contraindications, limitations, typical findings and possible side effects of breast tissue marker placement.
2. Store in a cool, dry place.
3. Do not use if the temperature indicator is black.
4. Use caution when handling the device to prevent premature deployment of the breast tissue marker.
5. Do not kink the flexible tube.
6. Maintain correct alignment of the yellow indicator (D in Figure 1) with the red arrow of the biopsy probe when dispensing pads.
7. Ensure that all pads are dispensed.
8. After use, this product may be a potential biohazard. Handle and dispose of in accordance with acceptable medical practice and applicable local, state, and federal laws and regulations.

POTENTIAL COMPLICATIONS

Complications may occur at any time during or after the procedure. Potential complications of breast tissue marker placement may include, but are not limited to: hematoma, hemorrhage, infection, adjacent tissue injury and pain.

EQUIPMENT REQUIRED


The following equipment is required to place the breast tissue marker:

- Appropriate imaging modality and accessories
- Appropriate biopsy probe
- Surgical gloves and drapes
- Local anesthetic
- Other equipment as necessary

HOW SUPPLIED

The SENOMARK[®] ULTRA Breast Tissue Marker is provided sterile and non-pyrogenic unless the package has been opened or damaged. **Sterilized by using irradiation. For single use only. Do not reuse. Do not resterilize.**


DIRECTIONS FOR USE (Refer to Figure 1)

1. Make certain that the collection of biopsy specimens has been completed. Ensure that the sample notch of the biopsy probe has been cleared of all tissue.
2. Press the "Marker or Anesthetic" button on the Control Module.

3. Remove the port cap from the sample container of the biopsy probe.
4. Insert the tip guide (B) into the sample container. Ensure that the tip guide is fully inserted. Refer to Figure 2.

Figure 2

5. Inspect the package to ensure that the package integrity has not been compromised. The product is sterile unless the seal is broken.
6. Using standard aseptic technique, remove the SENOMARK® ULTRA Breast Tissue Marker from the package and inspect for damage. REMOVE TIP PROTECTOR (A).
7. Insert the applicator through the tip guide into the biopsy probe.

Note: A small amount of resistance may be felt as the applicator handle (C) passes through the tip guide.

8. When the applicator contacts the proximal end of the probe, rotate the applicator, if necessary, and continue to advance the handle.
9. Advance the applicator to align the yellow indicator (D) with the red arrow of the biopsy probe. Refer to Figure 3.

Figure 3

Caution: To properly deploy pads, the yellow indicator must be fully seated against the probe and aligned with the red arrow.

10. Immediately deploy all pads by advancing the plunger (E).

Note: Avoid pulling the applicator handle during deployment which could cause deployment failure.

11. Rotate the biopsy probe 180° to position the sample notch away from the dispensed pads.
12. Remove the applicator from the biopsy probe.

Caution: If resistance is felt while removing the applicator, remove the entire probe/applicator assembly from the patient. Failure to do so may result in breakage of the applicator tip.

13. Dispose of the applicator properly.
14. Confirm final breast tissue marker position with imaging.

MR MR CONDITIONAL INFORMATION

Non-clinical testing has demonstrated that the breast tissue marker is MR Conditional. It can be scanned safely under the following conditions:

- Static magnetic field of 3-Tesla or less
- Spatial gradient field of 720-Gauss/cm or less
- Maximum whole-body-averaged specific absorption rate (SAR) of 3-W/kg for 15 minutes of scanning.

In non-clinical testing, the breast tissue marker produced a temperature rise of less than +0.4°C at a maximum whole body averaged specific absorption rate (SAR) of 3-W/kg of 15 minutes of MR scanning in a 3-Tesla, Signa MR system (Excite platform; G3.0-052B software, GE Healthcare, Milwaukee, WI).

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the breast tissue marker. Therefore, it may be necessary to optimize MR imaging parameters for the presence of this implant.

It is recommended that patients with a breast tissue marker register the MR conditions with the MedicAlert Foundation (www.medicalert.org).

WARRANTY

Bard Peripheral Vascular warrants to the first purchaser of this product that this product will be free from defects in materials and workmanship for a period of one year from the date of first purchase and liability under this limited product warranty will be limited to repair or replacement of the defective product, in Bard Peripheral Vascular's sole discretion or refunding your net price paid. Wear and tear from normal use or defects resulting from misuse of this product are not covered by this limited warranty.

TO THE EXTENT ALLOWABLE BY APPLICABLE LAW, THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL BARD PERIPHERAL VASCULAR BE LIABLE TO YOU FOR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM YOUR HANDLING OR USE OF THIS PRODUCT.

Some states/countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your state/country.

An issue or revision date and a revision number for these instructions are included for the user's information on the last page of this booklet. In the event 36 months have elapsed between this date and product use, the user should contact Bard Peripheral Vascular, Inc. to see if additional product information is available.

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Montované v Thajsku.

**SENO-MARK® ULTRA Breast Tissue Marker Ultrasound Enhanced Ribbon**

Marqueur de tissu mammaire SENOMARK® ULTRA pour aide échographique - ruban
 SENOMARK® ULTRA-Marker für Brustgewebe, durch Ultraschall verstärkte Schleife
 Nastro marcatore tissutale del tessuto mammario SENOMARK® ULTRA per eco-amplificazione
 Marcador de tejido mamario en forma de cinta SENOMARK® ULTRA para realce ecográfico
 SENOMARK® ULTRA borstweefselmarkeringshulpmiddel echografisch verbeterd lint
 Marcador de Tecido Mamário SENOMARK® ULTRA, Fita com Realce para Ecografia
 Υπερηχογραφικά ενισχυμένος ιχνηθέτης μαστικού ιστού SENOMARK® ULTRA τύπου ταινίας
 SENOMARK® ULTRA-bånd til forstærket ultralyd for brystvævsmarkør
 SENOMARK® ULTRA bröstvävnadsmarkör med ultraljudsförstärkande band
 SENOMARK® ULTRA -rintakudosmerkkilaitteen ultraäänivahvistusnauha
 SENOMARK® ULTRA-brystvevsmarkør, ultralydforsterket bånd
 Znacznik tkanki gruczołu sutkowego SENOMARK® ULTRA ze wstążką do znakowania ultrasonograficznego
 SENOMARK® ULTRA emlőszövet-jelölő ultrahangos jelzéssel ellátott szalag
 Páska pro zesílení ultrazvuku nástroje pro značení prsní tkáně SENOMARK® ULTRA
 SENOMARK® ULTRA Meme Doku İşaretleyicisi Ultrason Geliştirmeli Şerit
 Маркер места биопсии молочной железы SENOMARK® ULTRA — эхоконтрастная лента
 SENOMARK® ULTRA 乳房組織標記超音波增強型帶狀
 SENOMARK® ULTRA
 Ultrazvukový zosilňovač značkovača prsného tkaniva SENOMARK® ULTRA v tvare mašle

**SENO-MARK® ULTRA Breast Tissue Marker Ultrasound Enhanced Coil**

Marqueur de tissu mammaire SENOMARK® ULTRA pour aide échographique - bobine
 SENOMARK® ULTRA-Marker für Brustgewebe, durch Ultraschall verstärkte Wicklung
 Bobina marcatore tissutale del tessuto mammario SENOMARK® ULTRA per eco-amplificazione
 Marcador de tejido mamario en forma de bobina SENOMARK® ULTRA para realce ecográfico
 SENOMARK® ULTRA borstweefselmarkeringshulpmiddel echografisch verbeterde spoel
 Marcador de Tecido Mamário SENOMARK® ULTRA, Espiral com Realce para Ecografia
 Υπερηχογραφικά ενισχυμένος ιχνηθέτης μαστικού ιστού SENOMARK® ULTRA τύπου σπειράματος
 SENOMARK® ULTRA-spiral til forstærket ultralyd for brystvævsmarkør
 SENOMARK® ULTRA bröstvävnadsmarkör med ultraljudsförstärkande spiral
 SENOMARK® ULTRA -rintakudosmerkkilaitteen ultraäänivahvistuskierukka
 SENOMARK® ULTRA-brystvevsmarkør, ultralydforsterket spiral
 Znacznik tkanki gruczołu sutkowego SENOMARK® ULTRA ze zwojem do znakowania ultrasonograficznego
 SENOMARK® ULTRA emlőszövet-jelölő ultrahangos jelzéssel ellátott tekercs
 Cívka pro zesílení ultrazvuku nástroje pro značení prsní tkáně SENOMARK® ULTRA
 SENOMARK® ULTRA Meme Doku İşaretleyicisi Ultrason Geliştirmeli Bobin
 Маркер места биопсии молочной железы SENOMARK® ULTRA — эхоконтрастная катушка
 SENOMARK® ULTRA 乳房組織標記超音波增強型旋管狀
 SENOMARK® ULTRA
 Ultrazvukový zosilňovač značkovača prsného tkaniva SENOMARK® ULTRA v tvare cievky

**Contents**

Contenu
 Inhalt
 Contenuto
 Contenido
 Inhoud
 Conteúdo
 Περιεχόμενα
 Indhold
 Innehåll
 Sisältö
 Innhold
 Zawartość
 Tartalom
 Obsah
 İçindekiler
 Содержание
 裝物

Obsah

**Manufacturer**

Fabricant
 Hersteller
 Fabricante
 Fabricante
 Fabrikant
 Fabricante
 Κατασκευαστής
 Producent
 Tillverkare
 Valmistaja
 Fabrikant
 Producent
 Gyártó
 Výrobce
 Üretici
 Производитель
 製造商

Výrobca

25° C

Upper Limit of Temperature

Limite supérieure de température
 Obere Temperaturgrenze
 Limite superiore di temperatura
 Límite superior de temperatura
 Hoogste temperatuurlimiet
 Limite superior de temperatura
 Άνω όριο θερμοκρασίας
 Øvre temperaturgrænse
 Øvre temperaturgräns
 Lämpötilan yläraja
 Øvre temperaturgrense
 Górna wartość graniczna temperatury
 Felső hőmérséklet határ
 Horní hranice teploty
 Sıcaklık Üst Sınırı
 Верхний предел температуры
 溫度上限

Horný limit teploty

**Consult Instructions for Use**

Consulter le mode d'emploi
 Gebrauchsanweisung beachten
 Leggere le istruzioni per l'uso
 Consulte las Instrucciones de uso
 Raadpleeg de gebruiksaanwijzing
 Consultar as instruções de utilização
 Συμβουλευτείτε τις οδηγίες χρήσης
 Se brugsanvisningen
 Se bruksanvisning
 Lue käyttöohjeet
 Se bruksanvisning
 Zapoznać się z instrukcją obsługi
 Lásd a használati útmutatót
 Řiďte se návodem k použití
 Kullanma Talimatlarına Başvurun
 Обратитесь к инструкциям по применению
 請參考使用說明

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**Use By**

Date de péremption
 Verwendbar bis
 Data di scadenza
 Fecha de caducidad
 Uiterste gebruiksdatum
 Válido até
 Ημερομηνία λήξης
 Anvendes inden
 Användes före
 Käyttävänä ennen
 Bruk før
 Termin ważności
 Szavatosság lejár
 Datum expirace
 Son Kullanma Tarihi
 Срок годности
 有效期限

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**Gauge Size**

Diamètre
 Gauge-Maß
 Dimensioni gauge
 Tamaño de calibre
 Naaldidike
 Calibre
 Μέγεθος διαμετρήματος
 Gauge-størrelse
 Gauge-storlek
 G-koko
 Målerstørrelse
 Rozmiar igły
 Gauge méret
 Velikost měrky
 Kalibre Boyutu
 Калибр
 直徑尺寸

Menovitý rozmer

REF**Catalogue Number**

Número de catálogo
Katalognummer
Numero di catalogo
Número de referencia
Catalogusnummer
Número de Catálogo
Αριθμός καταλόγου
Katalognummer
Katalognummer
Luettelonumero
Katalognummer
Numer katalogowy
Katalógusszám
Katalogové číslo
Katalog Numarası
Номер по каталогу
目錄編號

Katalógové číslo

LOT**Lot Number**

Número de lot
Chargennummer
Numero di lotto
Número de lote
Lotnummer
Número do Lote
Αριθμός παρτίδας
Lotnummer
Lotnummer
Eränumero
Partinummer
Numer partii
Sorozatszám
Číslo šarže
Lot Numarası
Номер партии
批號
Lot
Číslo šarže

**Single Use**

À usage unique
Zur einmaligen Verwendung
Monouso
De un solo uso
Voor eenmalig gebruik
Não Reutilizar
Μίας χρήσης
Til engangsbrug
För engångsbruk
Kertakäyttöinen
Engangsbruk
Do jednorazowego użytku
Egyszer használatos
Pouze pro jednorázové použití
Tek Kullanılmıktır
Однократного применения
僅供一次性使用
1
Na jednorazové použitie

**Do Not Resterilize**

Ne pas résteriliser
Nicht erneut sterilisieren
Non risterilizzare
No volver a esterilizar
Niet opnieuw steriliseren
Não reesterilizar
Μην επαναποστεριώνετε
Må ikke reesteriliseres
Får ej reesteriliseras
Ei saa steriloida uudestaan
Må ikke steriliseres igjen.
Nie należy sterylizować ponownie
Újrasterilizálni tilos
Neresterilizujte
Tekrar Sterilize Etmeyin
Не подлежит повторной стерилизации
請勿重複消毒

Nesterilizujte opakovane

**Do Not Use if Temperature Indicator is Black**

Ne pas utiliser si l'indicateur thermique est noir.
Bei schwarz gefärbter Temperaturanzeige nicht verwenden
Non utilizzare se l'indicatore della temperatura è nero.
No usar si el punto de temperatura tiene color negro.
Niet gebruiken wanneer de temperatuursindicator zwart is.
Não utilize o indicador de temperatura estiver preto
Μη χρησιμοποιείτε εάν ο δείκτης θερμοκρασίας είναι μαύρος.
Må ikke anvendes, hvis temperaturpletten er sort
Använd inte produkten om temperaturindikatorn är svart.
Älä käytä, jos lämpötilan ilmaisin on musta
Skal ikke brukes hvis temperaturindikatoren er svart.
Nie używać, jeśli wskaźnik temperatury jest czarny
Ne használja, ha a hőmérséklet-kijelző fekete!
Produkt nepoužívejte, pokud je indikátor teploty černý
Sıcaklık Göstergesi Siyah Renklyse Kullanmayın
Запрещается использовать, если индикатор температуры имеет черный цвет.
如果溫度指示標籤是黑色的，請勿使用

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**Authorised Representative in the European Community**

Représentant autorisé au sein de la Communauté européenne
 Bevollmächtigter in der Europäischen Gemeinschaft
 Rappresentante autorizzato nella Comunità Europea
 Representante autorizado en la Comunidad Europea
 Gemachtigde binnen de Europese Gemeenschap
 Representante autorizado na Comunidade Europeia
 Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα
 Autoriseret repræsentant i Det Europæiske Fællesskab
 Auktoriserad representant inom EG
 Valtuutettu edustaja Euroopan yhteisössä
 Autorisert representant i EU
 Autoryzowany przedstawiciel na terenie Unii Europejskiej
 Meghatalmazott képviselő az Európai Közösségekben
 Autorizovaný zástupce v Evropské unii
 Avrupa Topluluğu Yetkili Temsilcisi
 Уполномоченный представитель в Европейском сообществе
 歐洲共同體的授權代表

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**Copyright ©2015 C. R. Bard, Inc. All rights reserved.**

Copyright ©2015 C. R. Bard, Inc. Tous droits réservés.
 Copyright ©2015 C. R. Bard, Inc. Alle Rechte vorbehalten.
 Copyright ©2015 C. R. Bard, Inc. Tutti i diritti riservati.
 Copyright ©2015 C. R. Bard, Inc. Todos los derechos reservados.
 Copyright ©2015 C.R. Bard, Inc. Alle rechten voorbehouden.
 Copyright ©2015 C. R. Bard, Inc. Todos os direitos reservados.
 Πνευματικά δικαιώματα ©2015. C. R. Bard, Inc. Με την επιφύλαξη
 αυτός δικαιώματος.
 Copyright ©2015 C. R. Bard, Inc. Alle rettigheder forbeholdt.
 Copyright ©2015 C. R. Bard, Inc. Med ensamrätt.
 Tekijänoikeus © 2015 C. R. Bard, Inc. Kaikki oikeudet pidätetään.
 Copyright ©2015 C. R. Bard, Inc. Med enerett.
 Copyright ©2015 C. R. Bard, Inc. Wszelkie prawa zastrzeżone.
 Copyright ©2015 C. R. Bard, Inc. Minden jog fenntartva.
 Copyright ©2015 C. R. Bard, Inc. Všechna práva vyhrazena.
 Telif hakkı ©2015 C. R. Bard, Inc. Tüm hakları saklıdır.
 Авторские права © 2015 г., C. R. Bard, Inc. Все права сохранены.
 Copyright ©2015 C. R. Bard, Inc. 保留所有權利。
 Copyright ©2015 C. R. Bard, Inc.
 Copyright ©2015 C. R. Bard, Inc. Všetky práva vyhrazené.

**Bard and SenoMark are trademarks and/or registered trademarks of C. R. Bard, Inc. or an affiliate. All other trademarks are the property of their respective owners.**

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**Authorised Representative
in the European Community**
Bard Limited
Forest House
Tilgate Forest Business Park
Brighton Road, Crawley
West Sussex
RH11 9BP, UK

BARID *Biopsy*

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Bard and SenoMark is C. R. Bard, Inc. or its affiliated organization's trademark and / or registered trademark. 所有其他商標為其所有人之所屬財產。
Bard SenoMark C. R. Bard, Inc.