SENO[®]MARK[®] ULTRA

SENO[®]MARK[®] ULTRA Breast Tissue Marker for Use with 9 or 12 Gauge ATEC[™] Biopsy Probe

INSTRUCTIONS FOR USE

CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

DESCRIPTION

The SENOMARK[®] ULTRA Breast Tissue Marker consists of a disposable applicator containing:

- 3 resorbable polyglycolic acid (PGA) pads that are essentially resorbed in approximately 12 weeks. The center pad contains a Titanium or BioDur[™] 108 wireform with an interwoven polyvinyl alcohol (PVA) polymer for permanent ultrasound enhancement. Please note the PVA polymer is not resorbed.
- 1 polyethylene glycol (PEG) push pellet.

The wireform is intended for long-term radiographic marking of the biopsy site. The PVA polymer is intended for long-term sonographic marking of the biopsy site.

The syringe-like applicator fits within the 9 or 12 gauge ATEC[™] Biopsy Probe to access the biopsy cavity. Pads are deployed from the applicator through the biopsy probe and into the biopsy cavity.

Product Code	Wireform Material	Wireform Shape	Probe Compatibility
SMAT9R	Titanium with PVA	Ribbon 	9G
SMAT9C	BioDur [™] 108 with PVA	Coil 	9G
SMAT12R	Titanium with PVA	Ribbon 	12G
SMAT12C	BioDur [™] 108 with PVA	Coil 	12G

INDICATIONS FOR USE

The SENOMARK[®] ULTRA Breast Tissue Marker is intended to radiographically and sonographically mark breast tissue during a percutaneous breast biopsy procedure.

CONTRAINDICATIONS

Patients with a known hypersensitivity to the materials listed in the device description may suffer an allergic reaction to this implant.

WARNINGS

1. Use caution when placing near a breast implant to avoid puncture of the implant capsule.
2. Avoid the use of excessive force during removal of the applicator to prevent breakage of the applicator tip.
3. Do not use in the presence of infection.
4. This device has been designed for single use only. Reusing this medical device bears the risk of cross-patient contamination as medical devices – particularly those with long and small lumina, joints, and/or crevices between components – are difficult or impossible to clean once body fluids or tissues with potential pyrogenic or microbial contamination have had contact with the medical device for an indeterminable period of time. The residue of biological material can promote the contamination of the device with pyrogens or microorganisms which may lead to infectious complications.
5. Do not resterilize. After reesterilization, the sterility of the product is not guaranteed because of an indeterminable degree of potential pyrogenic or microbial contamination which may lead to infectious complications. Cleaning, reprocessing and/or reesterilization of the present medical device increases the probability that the device will malfunction due to potential adverse effects on components that are influenced by thermal and/or mechanical changes.

PRECAUTIONS

1. This product should only be used by a physician who is completely familiar with the indications, contraindications, limitations, typical findings and possible side effects of breast tissue marker placement.
2. Store in a cool, dry place.
3. Do not use if the temperature indicator is black.
4. Use caution when handling the device to prevent premature deployment of the breast tissue marker.
5. Do not kink the flexible tube.
6. Maintain correct alignment of the alignment indicator (D in Figure 1) with the sample notch when dispensing pads.
7. Ensure that all pads are dispensed.
8. After use, this product may be a potential biohazard. Handle and dispose of in accordance with acceptable medical practice and applicable local, state, and federal laws and regulations.

POTENTIAL COMPLICATIONS

Complications may occur at any time during or after the procedure. Potential complications of breast tissue marker placement may include, but are not limited to: hematoma, hemorrhage, infection, adjacent tissue injury and pain.

EQUIPMENT REQUIRED


The following equipment is required to place the breast tissue marker:

- Appropriate imaging modality and accessories
- Appropriate biopsy probe
- Surgical gloves and drapes
- Local anesthetic
- Other equipment as necessary

HOW SUPPLIED

The SENOMARK[®] ULTRA Breast Tissue Marker is provided sterile and non-pyrogenic unless the package has been opened or damaged. **Sterilized by using irradiation. For single use only. Do not reuse. Do not resterilize.**

DIRECTIONS FOR USE (Refer to Figure 1)

1. Make certain that the collection of biopsy specimens has been completed. Ensure that the sample notch of the biopsy probe has been cleared of all tissue.
2. Inspect the package to ensure that the package integrity has not been compromised. The product is sterile unless the seal is broken.
3. Remove the sample collection chamber from the back of the ATEC[™] Biopsy Probe. Use the provided cap (B) to seal the collection chamber after removing the chamber from the biopsy probe.
4. Remove the guide seal (C) from the package and lock into the back of the biopsy probe.
5. Using standard aseptic technique, remove the SENOMARK[®] ULTRA Breast Tissue Marker from the package and inspect for damage. REMOVE TIP PROTECTOR (A).

6. Insert the applicator into the biopsy probe through the guide seal. Advance the applicator until the applicator contacts the distal end of the biopsy probe.
7. Align the alignment indicator (D) on the applicator handle (E) with the sample notch in the biopsy probe.
8. Immediately deploy all pads by advancing the plunger (F). The pads will be fully dispensed when the plunger contacts the handle.

Note: Avoid pulling the applicator handle during deployment which could cause deployment failure.

9. Rotate the biopsy probe 180° to position the sample notch away from the dispensed pads.
10. Remove the applicator from the biopsy probe. Avoid the use of force.

Caution: If resistance is felt while removing the applicator, remove the entire probe/applicator assembly from the patient. Failure to do so may result in breakage of the applicator tip.

11. Close the sample notch and remove the biopsy probe, per manufacturer's instructions.
12. Dispose of the applicator properly.
13. Confirm final breast tissue marker position with imaging.

MR CONDITIONAL INFORMATION

Non-clinical testing has demonstrated that the breast tissue marker is MR Conditional. It can be scanned safely under the following conditions:

- Static magnetic field of 3-Tesla or less
- Spatial gradient field of 720-Gauss/cm or less
- Maximum whole-body-averaged specific absorption rate (SAR) of 3-W/kg for 15 minutes of scanning.

In non-clinical testing, the breast tissue marker produced a temperature rise of less than +0.4°C at a maximum whole body averaged specific absorption rate (SAR) of 3-W/kg of 15 minutes of MR scanning in a 3-Tesla, Signa MR system (Excite platform; G3.0-052B software, GE Healthcare, Milwaukee, WI).

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the breast tissue marker. Therefore, it may be necessary to optimize MR imaging parameters for the presence of this implant.

It is recommended that patients with a breast tissue marker register the MR conditions with the MedicAlert Foundation (www.medicalert.org).

WARRANTY

Bard Peripheral Vascular warrants to the first purchaser of this product that this product will be free from defects in materials and workmanship for a period of one year from the date of first purchase and liability under this limited product warranty will be limited to repair or replacement of the defective product, in Bard Peripheral Vascular's sole discretion or refunding your net price paid. Wear and tear from normal use or defects resulting from misuse of this product are not covered by this limited warranty.

TO THE EXTENT ALLOWABLE BY APPLICABLE LAW, THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL BARD PERIPHERAL VASCULAR BE LIABLE TO YOU FOR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM YOUR HANDLING OR USE OF THIS PRODUCT.

Some states/countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your state/country.

An issue or revision date and a revision number for these instructions are included for the user's information on the last page of this booklet. In the event 36 months have elapsed between this date and product use, the user should contact Bard Peripheral Vascular, Inc. to see if additional product information is available.

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Bard и SenoMark are C. R. Bard, Inc. or its affiliated company's trademarks and/or registered trademarks. All other trademarks are the property of their respective owners.

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**Authorised Representative in the European Community**

Représentant autorisé au sein de la Communauté européenne
 Bevollmächtigter in der Europäischen Gemeinschaft
 Rappresentante autorizzato nella Comunità Europea
 Representante autorizado en la Comunidad Europea
 Gemachtigde binnen de Europese Gemeenschap
 Representante autorizado na Comunidade Europeia
 Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα
 Autoriseret repræsentant i Det Europæiske Fællesskab
 Auktoriserad representant inom EG
 Valtuutettu edustaja Euroopan yhteisössä
 Autoriseret representant i EU
 Autoryzowany przedstawiciel na terenie Unii Europejskiej
 Meghatalmazott képviselő az Európai Közösségben
 Autorizovaný zástupce v Evropské unii
 Avrupa Topluluğu Yetkili Temsilcisi
 Уполномоченный представитель в Европейском сообществе
 歐洲共同體的授權代表

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**Copyright ©2015 C. R. Bard, Inc. All rights reserved.**

Copyright ©2015 C. R. Bard, Inc. Tous droits réservés.
 Copyright ©2015 C. R. Bard, Inc. Alle Rechte vorbehalten.
 Copyright ©2015 C. R. Bard, Inc. Tutti i diritti riservati.
 Copyright ©2015 C. R. Bard, Inc. Todos los derechos reservados.
 Copyright ©2015 C.R. Bard, Inc. Alle rechten voorbehouden.
 Copyright ©2015 C. R. Bard, Inc. Todos os direitos reservados.
 Πνευματικά δικαιώματα ©2015. C. R. Bard, Inc. Με την επιφύλαξη παντός δικαιώματος.
 Copyright ©2015 C. R. Bard, Inc. Alle rettigheder forbeholdt.
 Copyright ©2015 C. R. Bard, Inc. Med ensamrätt.
 Tekijänoikeus © 2015 C. R. Bard, Inc. Kaikki oikeudet pidätetään.
 Copyright ©2015 C. R. Bard, Inc. Med enerett.
 Copyright ©2015 C. R. Bard, Inc. Wszelkie prawa zastrzeżone.
 Copyright ©2015 C. R. Bard, Inc. Minden jog fenntartva.
 Copyright ©2015 C. R. Bard, Inc. Všechna práva vyhrazena.
 Telif hakkı ©2015 C. R. Bard, Inc. Tüm hakları saklıdır.
 Авторские права © 2015 r., C. R. Bard, Inc. Все права сохранены.
 Copyright ©2015 C. R. Bard, Inc. 保留所有權利。
 Copyright ©2015 C. R. Bard, Inc.
 Copyright ©2015 C. R. Bard, Inc. Všetky práva vyhrazené.

Manufacturer:
Bard Peripheral Vascular, Inc.
1625 West 3rd Street
Tempe, AZ 85281
USA
TEL: 1-480-894-9515
1-800-321-4254
FAX: 1-480-966-7062
1-800-440-5378
www.bardbiopsy.com

**Authorised Representative
in the European Community**
Bard Limited
Forest House
Tilgate Forest Business Park
Brighton Road, Crawley
West Sussex
RH11 9BP, UK

BARID *Biopsy*

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Bard and SenoMark are trademarks and/or registered trademarks of C. R. Bard, Inc. or an affiliate. All other trademarks are the property of their respective owners.