SENO^{MARK}
SENO^{MARK} Biopsy Site Marker
For Use With ENCOR[®] Biopsy Probes

ENGLISH

- REF** SMTEC10G
 SMEC10GSS
 SMEC7GSS
 SMEC12GSS

INSTRUCTIONS FOR USE

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.


Figure 1 - Applicator

DESCRIPTION

The SENOMARK[®] Biopsy Site Marker consists of a disposable applicator containing:

- 3 resorbable polyglycolic acid (PGA) pads. The center pad contains a titanium or stainless steel wireform.
- 1 polyethylene glycol (PEG) push pad at the proximal end.

Product Code	Wireform Material	Wireform Shape	Probe Compatibility
SMTEC10G	Titanium		10 G ENCOR [®]
SMEC10GSS	316L stainless steel		10 G ENCOR [®]
SMEC7GSS	316L stainless steel		7 G ENCOR [®]
SMEC12GSS	316L stainless steel		12 G ENCOR [®]

The wireform is intended for long-term radiographic marking of the biopsy site. The pads are visible via ultrasound for approximately 3 weeks and are essentially resorbed in approximately 12 weeks.

The syringe-like applicator fits within the ENCOR[®] probe to access the biopsy cavity. Pads are deployed from the Applicator through the biopsy device into the biopsy cavity.

INTENDED USE

The SENOMARK[®] Biopsy Site Marker is intended to radiographically and sonographically mark breast tissue during a percutaneous breast biopsy procedure.

CONTRAINDICATIONS

None known

WARNINGS

- This device has been designed for single use only. Reusing this medical device bears the risk of cross-patient contamination as medical devices – particularly those with long and small lumina, joints, and/or crevices between components – are difficult or impossible to clean once body fluids or tissues with potential pyrogenic or microbial contamination have had contact with the medical device for an indeterminable period of time. The residue of biological material can promote the contamination of the device with pyrogens or microorganisms which may lead to infectious complications.
- Do not resterilize. After resterilization, the sterility of the product is not guaranteed because of an indeterminable degree of potential pyrogenic or microbial contamination which may lead to infectious complications. Cleaning, reprocessing and/or resterilization of the present medical device increases the probability that the device will malfunction due to potential adverse effects on components that are influenced by thermal and/or mechanical changes.
- Patients with a known hypersensitivity to the materials listed in the device description may suffer an allergic reaction to this implant.
- The SENOMARK[®] is not recommended for use in patients with breast implants.
- Do not use in the presence of infection.
- Avoid the use of excessive force during removal of the Applicator to prevent breakage of the Applicator tip.

PRECAUTIONS

- The device should only be used by physicians trained in percutaneous biopsy procedures.
- Do not use this product if the sterile barrier has been previously opened or if the package is damaged.
- This product is provided sterile, and is intended for single use only. Do not resterilize.
- Do not kink the flexible tube.
- Maintain correct alignment of the yellow indicator key with the red arrow of the biopsy probe when dispensing pads (Figure 1).
- Ensure that all pads are dispensed.
- Published studies with comparably sized 316L stainless steel biopsy markers have shown no marker movement and insignificant marker heating when tested in an MRI System with a 1.5T field strength.
- Store at temperatures below 25°C (77°F).
- Do not use if temperature dot is black.

COMPLICATIONS

Potential complications that may be associated with the use of the SENOMARK[®] biopsy site marker are similar to those associated with the use of other biopsy marking devices.

DIRECTIONS FOR USE

Refer to Figure 1 for Applicator diagram.

1. Perform a "dry tap" to ensure that the sample notch is clear of tissue.
2. Using standard aseptic technique, remove the Applicator from the package and inspect for damage. Remove Tip Protector (A).
3. Press the "Marker or Anesthetic" button on the Control Module.
4. Remove the Port Cap from of the Sample Chamber of the probe.
5. Insert Tip Guide (B) into Sample Chamber. Ensure that the Tip Guide is fully inserted. Refer to Figure 2.


Figure 2: Insert Tip Guide into Sample Chamber

6. Remove the Tip Protector from the Applicator.
7. Insert the Applicator through the Tip Guide (B) into the biopsy probe.
Note: A small amount of resistance may be felt as the Applicator Handle (C) passes through the Tip Guide.
8. When the Applicator contacts the proximal end of the probe, rotate the Applicator, if necessary, and continue to advance the handle.
9. Advance the Applicator to align the Yellow Indicator with the red arrow of the biopsy probe. Refer to Figure 3.


Figure 3: Applicator alignment in probe

Caution: To properly deploy pads, the Applicator Yellow Indicator must be fully seated against the probe and aligned with the red arrow. See Figure 4.


Figure 4: Fully seated Applicator in probe

10. Immediately deploy all pads by advancing the Plunger (E).
 - Note:** Avoid pulling the Applicator Handle during deployment which could cause deployment failure.
11. Rotate the biopsy device 180° to position the Sample Chamber away from the dispensed pads.
12. Remove the Applicator from the biopsy device.
 - Caution:** If resistance is felt while removing the Applicator, remove the entire probe/Applicator assembly from the patient. Failure to do so may result in breakage of the Applicator Tip.
13. Dispose of the Applicator properly.
14. Confirm final marker position with imaging.

WARRANTY

Bard Peripheral Vascular warrants to the first purchaser of this product that this product will be free from defects in materials and workmanship for a period of one year from the date of first purchase and liability under this limited product warranty will be limited to repair or replacement of the defective product, in Bard Peripheral Vascular's sole discretion or refunding your net price paid. Wear and tear from normal use or defects resulting from misuse of this product are not covered by this limited warranty.

TO THE EXTENT ALLOWABLE BY APPLICABLE LAW, THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL BARD PERIPHERAL VASCULAR BE LIABLE TO YOU FOR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM YOUR HANDLING OR USE OF THIS PRODUCT.

Some states/countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your state/country.

An issue or revision date and a revision number for these instructions are included for the user's information on the last page of this booklet. In the event 36 months have elapsed between this date and product use, the user should contact Bard Peripheral Vascular, Inc. to see if additional product information is available.

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Copyright © 2015 C.R. Bard, Inc. All rights reserved.
Copyright © 2015 C.R. Bard, Inc. Tous droits réservés.
Copyright © 2015 C.R. Bard, Inc. Alle Rechte vorbehalten.
Copyright © 2015 C.R. Bard, Inc. Tutti i diritti riservati.
Copyright © 2015 C.R. Bard, Inc. Todos los derechos reservados.
Copyright © 2015 C.R. Bard, Inc. Alle rechten voorbehouden.