SENO[®]MARK[®]
SENO[®]MARK[®] Biopsy Site Marker
For Use with a 9 or 12 Gauge ATEC[™]

ENGLISH

REF SMRSU9GT
 SMRSU12GT

INSTRUCTIONS FOR USE

Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

DESCRIPTION

The SENOMARK[®] Biopsy Site Marker consists of a disposable applicator containing:
 3 resorbable polyglycolic acid (PGA) pads, the center pad contains a titanium wireform.
 1 polyethylene glycol (PEG) push pellet.

Product Code	Wireform Material	Wireform Shape	Probe Compatibility
SMRSU9GT	Titanium	X	9G ATEC [™]
SMRSU12GT	Titanium	S	12G ATEC [™]

The wireform is intended for long-term radiographic marking of the biopsy site. The pads are visible via ultrasound for approximately 3 weeks and are essentially resorbed in approximately 12 weeks.

The syringe-like applicator fits within the 9 or 12 gauge ATEC[™] to access the biopsy cavity¹. Pads are deployed from the applicator through the biopsy device and into the biopsy cavity.

INTENDED USE

The SENOMARK[®] is intended to radiographically and sonographically mark breast tissue during a percutaneous breast biopsy procedure using a 9 or 12 gauge Suros ATEC[™] biopsy device.

CONTRAINDICATIONS

This device is not intended for use except as indicated.

WARNINGS

- This device has been designed for single use only. Reusing this medical device bears the risk of cross-patient contamination as medical devices – particularly those with long and small lumina, joints, and/or crevices between components – are difficult or impossible to clean once body fluids or tissues with potential pyrogenic or microbial contamination have had contact with the medical device for an indeterminable period of time. The residue of biological material can promote the contamination of the device with pyrogens or microorganisms which may lead to infectious complications.
- Do not resterilize. After resterilization, the sterility of the product is not guaranteed because of an indeterminable degree of potential pyrogenic or microbial contamination which may lead to infectious complications. Cleaning, reprocessing and/or resterilization of the present medical device increases the probability that the device will malfunction due to potential adverse effects on components that are influenced by thermal and/or mechanical changes.
- Patients with a known hypersensitivity to the materials listed in the device description may suffer an allergic reaction to this implant.
- The SENOMARK[®] is not recommended for use in patients with breast implants.
- Do not use in the presence of infection.

PRECAUTIONS

- The device should only be used by physicians trained in percutaneous biopsy procedures.
- Do not use this product if the sterile barrier has been previously opened or if the package is damaged.
- This product is provided sterile, and is intended for single use only. **DO NOT RESTERILIZE.**
- Maintain correct orientation of the Alignment Indicator (A) with the sample notch when dispensing pads.
- Ensure that all pads are dispensed.
- Store at temperatures below 25° C (77° F).

COMPLICATIONS


Potential complications (e.g. infection) that may be associated with the use of the SENOMARK[®] are the same as those associated with the use of other biopsy marking devices.

HOW SUPPLIED

The SENOMARK[®] is provided sterile, and is intended for single use only.

DIRECTIONS FOR USE (Refer to Figure)

1. Make certain that the collection of biopsy specimens has been completed. Ensure that the sample notch has been cleared of all tissue.
2. Inspect the package to ensure that the package integrity has not been compromised. The product is sterile unless the seal is broken.
3. Remove the sample collection chamber from the back of the ATEC[™] Handpiece. Use the provided Cap (B) to seal the collection chamber after removing the chamber from the biopsy device.
4. Remove the Guide Seal (C) from the package and lock into the back of the ATEC[™] Handpiece.
5. Using standard aseptic technique, remove the SENOMARK[®] device from the package and check for damage. REMOVE TIP PROTECTOR (D).
6. Place the SENOMARK[®] applicator into the ATEC[™] biopsy device through the Guide Seal. Advance the applicator until the applicator contacts the distal end of the ATEC[™] cannula.
7. Align the SENOMARK[®] Alignment Indicator (A) on the Applicator Handle (E) with the sample notch in the biopsy device cannula.
8. Immediately deploy the SENOMARK[®] pads by advancing the Plunger (F) while maintaining forward pressure on the Applicator Handle (E), and correct alignment of the applicator with the sample notch. The pads will be fully dispensed when the Plunger contacts the Handle.
9. Rotate the biopsy device 180 degrees to position the sample notch away from the dispensed pellets.
10. Remove the SENOMARK[®] applicator from the biopsy device cannula. Avoid the use of force.
11. Close the sample notch and remove the biopsy device, per manufacturer's instructions.
12. Dispose of the applicator properly.
13. Confirm final marker position with imaging.

SENO[®]MARK[®] Biopsy Marker

WARRANTY

Bard Peripheral Vascular warrants to the first purchaser of this product that this product will be free from defects in materials and workmanship for a period of one year from the date of first purchase and liability under this limited product warranty will be limited to repair or replacement of the defective product, in Bard Peripheral Vascular's sole discretion or refunding your net price paid. Wear and tear from normal use or defects resulting from misuse of this product are not covered by this limited warranty.

TO THE EXTENT ALLOWABLE BY APPLICABLE LAW, THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL BARD PERIPHERAL VASCULAR BE LIABLE TO YOU FOR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM YOUR HANDLING OR USE OF THIS PRODUCT.

Some states/countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your state/country.

An issue or revision date and a revision number for these instructions are included for the user's information on the last page of this booklet. In the event 36 months have elapsed between this date and product use, the user should contact Bard Peripheral Vascular, Inc. to see if additional product information is available.

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
**Do Not Use if Temperature Indicator is Black**

Ne pas utiliser si l'indicateur thermique est noir.
 Bei schwarz gefärbter Temperaturanzeige nicht verwenden
 Non utilizzare se l'indicatore della temperatura è nero.
 No usar si el punto de temperatura tiene color negro.
 Niet gebruiken wanneer de temperatuursindicator zwart is.
 Não utilize o indicador de temperatura estiver preto
 Μη χρησιμοποιείτε εάν ο δείκτης θερμοκρασίας είναι μαύρος.
 Må ikke anvendes, hvis temperaturpletten er sort
 Använd inte produkten om temperaturindikatorn är svart.
 Älä käytä, jos lämpötilan ilmaisin on musta
 Skäl ikke brukes hvis temperaturindikatoren er svart.
 Nie używać, jeśli wskaźnik temperatury jest czarny
 Ne használja, ha a hőmérséklet-kijelző feketé!
 Produkt nepoužívejte, pokud je indikátor teploty černý
 Sıcaklık Göstergesi Siyah Renkliye Kullanmayın
 Запрещается использовать, если индикатор температуры имеет черный цвет.
 如果溫度指示標籤是黑色的，請勿使用

†

**Gauge Size**

Diamètre
 Gauge-Maß
 Dimensiuni gauge
 Tamaño de calibre
 Naaldidikte
 Calibre
 Μέγεθος διαμετρήματος
 Gauge-størrelse
 Gauge-storlek
 G-koko
 Målerstørrelse
 Rozmiar igły
 Gauge méret
 Velikost měřky
 Kalibre Boyutu
 Калибр
 直徑尺寸

**Upper Limit of Temperature**

Limite supérieure de température
 Obere Temperaturgrenze
 Temperatura massima
 Límite superior de temperatura
 Bovengrens temperatuur
 Limite superior de temperatura
 Ανώτατο όριο θερμοκρασίας
 Øvre temperaturgrense
 Högsta temperatur
 Lämpötilan yläraja
 Øvre temperaturgrense
 Maksymalna dopuszczalna temperatura
 A hőmérséklet felső határértéke
 Horní mez teploty
 Üst Sıcaklık Limiti
 Верхний предел температуры
 溫度上限

**Not Made with Natural Rubber Latex**

Fabriqu e sans latex naturel
 Nicht aus Naturgummilatem hergestellt.
 Non prodotto con lattice di gomma naturale
 Este producto no se fabrica con latex de caucho natural
 Niet gemaakt met latex van natuurlijk rubber
 N o cont em latex de borracha natural
 Δεν κατασκευάζεται από φυσικό ελαστικό λάτεξ
 Ikke fremstillet med naturlig gummilatem
 Ej tillverkad med naturgummilatem
 Valmistuksessa ei ole k ytetty luonnonkumilatemksia
 Dette produktet er ikke produsert med naturlig gummilatemks
 Produkt nie jest wytwarzany z lateksu kauczuku naturalnego
 Nem tartalmaz természetes gumilatemet.
 Nenı v vyroben z p ırodnıho pry ov eho latexu
 Dođal Kau uk Lateksten  retimemiştir
 He содержит натуральный каучуковый латекс
 製造未採用天然膠乳

**Sterilized Using Irradiation**

St rilis e par irradiation
 Durch Bestrahlung sterilisiert
 Sterilizzato mediante irradiazione
 Esterilizado mediante irradiaci n
 Gesteriliseerd met behulp van straling
 Esterilizado por Irradia o
 Αποστειρωμένο με χρήση ακτινοβολίας
 Steriliseret ved brug af bestr ling
 Steriliserad med str lning
 Steriloitu s teilytt m ll 
 Steriliseret ved str lning
 Produkt sterylizowany przez napromienianie
 Sug rzs ssal steriliz lt
 Sterilizov no oz a enim
 Radyoaktif Işınla Sterilize Edilmiştir
 Стерилизовано радиацией
 已使用放射線消毒

**Do Not Use If Package Is Damaged**

Ne pas utiliser si l'emballage est endommag e
 Nicht verwenden, wenn die Verpackung besch digt ist
 Non utilizzare se la confezione   danneggiata
 No usar si el envase est  da ado
 Niet gebruiken wanneer de verpakking beschadigd is
 N o utilizar se a embalagem estiver danificada
 Μη χρησιμοποιείτε το πακ on εάν η συσκευασία του  χει υποστεί ζημι 
 M  ikke bruges, hvis emballagen er beskadiget
 Anv nd inte produkten om f rpackningen skadats
 Ei saa k ytt t , jos pakkaus on vaurioitunut
 Sk l ikke brukes hvis pakningen er skadet
 Nie stosowa , jeřli opakowanie jest uszkodzone
 Ne haszn lja fel, ha a csomagol s s r lt
 Nepou ıvejte, pokud je obal pořkozen
 Paket Hasarlıya Kullanmayınız
 He использовать, если упаковка повреждена
 如果包裝已受損，請勿使用

**Breast Tissue Marker**

Marqueur de tissu mammaire
 Brustgewebemarker
 Marcatore tessutale mammario
 Marcador de tejido mamario
 Mammaweefselmarkering
 Marcador de tecido mam rio
 Δείκτης μαστικής ιστού
 Brystv vnsmark r
 Br stv vnadsmark r
 Rintakudosmarkkeri
 Brystvevmark r
 Marker tkanki piersiowej
 Eml sz vet marker
 Zna kova  prsnı tk n 
 Meme Dokusu İřaretleyicisi
 Биологический маркер молочной железы
 乳房組織標記

**Authorised Representative in the European Community**

Repr sentant autoris  au sein de la Communaut  europ enne
 Berechtigter Vertreter in der Europ ischen Gemeinschaft
 Rappresentante autorizzato nella Comunit  Europea
 Representante autorizado en la Comunidad Europea
 Gemachtigde vertegenwoordiger binnen de Europese Gemeenschap
 Representante autorizado na Comunidade Europeia
 Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα
 Autoriseret repr sentant i Det Europ iske F llesskab
 Auktoriserad representant i EG
 Valtuutettu edustaja Euroopan yhteis ss 
 Autoriseret representant i EU
 Autoryzowany przedstawiciel na terenie Unii Europejskiej
 Bejegyzett k pviselet az Eur pai K z ss gben
 Autorizovaný z stupce pro Evropskou unii
 Avrurpa Topluļuğu Yetkili Temsilcisi
 Уполномоченный представитель в Европейском сообществе
 歐洲共同體的授權代表

Copyright © 2015 C.R. Bard, Inc. All rights reserved.

Copyright © 2015 C.R. Bard, Inc. Tous droits réservés.
Copyright © 2015 C.R. Bard, Inc. Alle Rechte vorbehalten.
Copyright © 2015 C.R. Bard, Inc. Tutti i diritti riservati.
Copyright © 2015 C.R. Bard, Inc. Todos los derechos reservados.
Copyright © 2015 C.R. Bard, Inc. Alle rechten voorbehouden.
Copyright © 2015 C.R. Bard, Inc. Todos os direitos reservados.
Πνευματικά δικαιώματα © 2015 C.R. Bard, Inc. Με την επιφύλαξη παντός δικαιώματος.
Copyright © 2015 C.R. Bard, Inc. Alle rettigheder forbeholdt.
Copyright © 2015 C.R. Bard, Inc. Eftertryck förbjudes.
Copyright © 2015 C.R. Bard, Inc. Kaikki oikeudet pidätetään.
Copyright © 2015 C.R. Bard, Inc. Med enerrett.
Copyright © 2015 C.R. Bard, Inc. Wszelkie prawa zastrzeżone.
Copyright © 2015 C.R. Bard, Inc. Minden jog fenntartva.
Copyright © 2015 C.R. Bard, Inc. Všechna práva vyhrazena.
Telif Hakki © 2015 C.R. Bard, Inc. Tüm hakları saklıdır.
Авторские права © 2015 Компания С. R. Bard, Inc. Все права защищены.
Copyright © 2015 C.R. Bard, Inc. 保留所有權利。
© 2015 C.R. Bard, Inc.

Bard and SenoMark are trademarks and/or registered trademarks of C. R. Bard, Inc. or an affiliate. All other trademarks are the property of their respective owners.

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 **Manufacturer:**
SenoRx, Inc.
1625 West 3rd Street
Tempe, AZ 85281
USA

TEL: 1-480-894-9515
1-800-321-4254
FAX: 1-480-966-7062
1-800-440-5376
www.bardbiopsy.com

0086

**Authorized Representative in
the European Community**
Bard Limited
Forest House
Tilgate Forest Business Park
Brighton Road, Crawley
West Sussex
RH11 9BP, UK

BARB *Biopsy*

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Bard SenoMark C. R. Bard, Inc.