SENOMARK®

SENOMARK® Biopsy Site Marker


For Use With Eviva™ 9 Ga Biopsy Probes

ENGLISH

REF SMSE9GT

INSTRUCTIONS FOR USE

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

DESCRIPTION

The SENOMARK® Biopsy Site Marker consists of a disposable applicator containing: 3 resorbable polyglycolic acid (PGA) pads. The Center pad contains a titanium wireform. 1 polyethylene glycol (PEG) push pellet proximal of the PGA pads.

Product Code	Wireform Material	Wireform Shape	Probe Compatibility
SMSE9GT	Titanium	X	9G Eviva™

The wireform is intended for long term radiographic marking of the biopsy site. The pads are visible via ultrasound for approximately 3 weeks and are essentially resorbed in approximately 12 weeks.

The syringe-like applicator fits within the Eviva™ cannula to access the biopsy cavity pads are deployed from the applicator through the biopsy device into the biopsy cavity.

INTENDED USE

The SENOMARK® is intended to radiographically and ultrasonically mark breast tissue during a percutaneous biopsy procedure.

CONTRAINDICATIONS

This device is not intended for use except as indicated.

WARNINGS

- This device has been designed for single use only. Reusing this medical device bears the risk of cross-patient contamination as medical devices – particularly those with long and small lumina, joints, and/or crevices between components – are difficult or impossible to clean once body fluids or tissues with potential pyrogenic or microbial contamination have had contact with the medical device for an indeterminable period of time. The residue of biological material can promote the contamination of the device with pyrogens or microorganisms which may lead to infectious complications.
- Do not resterilize. After resterilization, the sterility of the product is not guaranteed because of an indeterminable degree of potential pyrogenic or microbial contamination which may lead to infectious complications. Cleaning,

reprocessing and/or resterilization of the present medical device increases the probability that the device will malfunction due to potential adverse effects on components that are influenced by thermal and/or mechanical changes.

- Patients with a known hypersensitivity to the materials listed in the device description may suffer an allergic reaction to this implant.
- The SENOMARK® is not recommended for use in patients with breast implants.
- Do not use in the presence of infection.

PRECAUTIONS

- The device should only be used by physicians trained in percutaneous biopsy procedures.
- Do not use this product if the sterile barrier has been breached.
- This product is provided sterile, and is intended for single use only.
- DO NOT RESTERILIZE.**
- Special care should be taken not to kink the flexible tube.
- Published studies with comparably sized 316L stainless steel biopsy markers have shown no marker movement and insignificant marker heating when tested in an MRI System with a 1.5T field strength.
- Store at temperatures below 25°C (77°F).

COMPLICATIONS

Potential complications (e.g. infection) that may be associated with the use of the SENOMARK® are similar to those associated with the use of other biopsy marking devices.

HOW SUPPLIED

The SENOMARK® Eviva™ is provided sterile, and is intended for single use only.

DIRECTIONS FOR USE (Refer to Figure 1)

- Make sure that biopsy specimen collection has been completed.
- While holding the cannula in place, remove the biopsy device from the cannula (refer to Eviva™ Biopsy Probe instructions for use).
- Inspect the package to ensure that the package integrity has not been compromised.
- Using standard aseptic technique, remove the SENOMARK® device from the package and check for damage. Remove tip cover.
- While holding the biopsy device cannula in place, insert the applicator into the cannula and advance it until the applicator stop (B) contacts the cannula.
- Immediately deploy the pads by advancing the applicator plunger (A) while maintaining forward pressure on the handle (C). The pads will be fully dispensed when the plunger contacts the finger tabs of the handle (C).
- Remove the applicator and biopsy device cannula. Dispose of property.
- Confirm final marker position with imaging.

WARRANTY

Bard Peripheral Vascular warrants to the first purchaser of this product that this product will be free from defects in materials and workmanship for a period of one year from the date of first purchase and liability under this limited product warranty will be limited to repair or replacement of the defective product, in Bard Peripheral Vascular's sole discretion or refunding your net price paid. Wear and tear from normal use or defects resulting from misuse of this product are not covered by this limited warranty.

TO THE EXTENT ALLOWABLE BY APPLICABLE LAW, THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL BARD PERIPHERAL VASCULAR BE LIABLE TO YOU FOR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM YOUR HANDLING OR USE OF THIS PRODUCT.

Some states/countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your state/country.

An issue or revision date and a revision number for these instructions are included for the user's information on the last page of this booklet. In the event 36 months have elapsed between this date and product use, the user should contact Bard Peripheral Vascular, Inc. to see if additional product information is available.

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Copyright © 2015 C.R. Bard, Inc. All rights reserved.
Copyright © 2015 C.R. Bard, Inc. Tous droits réservés.
Copyright © 2015 C.R. Bard, Inc. Alle Rechte vorbehalten.
Copyright © 2015 C.R. Bard, Inc. Tutti i diritti riservati.
Copyright © 2015 C.R. Bard, Inc. Todos los derechos reservados.
Copyright © 2015 C.R. Bard, Inc. Alle rechten voorbehouden.