# SENO<sup>®</sup>MARK<sup>®</sup>

## SENO<sup>®</sup>MARK<sup>®</sup> Biopsy Site Marker

For use with 11 GAUGE Mammotome™

**ENGLISH**

**REF** SMTMT11G

### INSTRUCTIONS FOR USE

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

### DESCRIPTION

The SENOMARK<sup>®</sup> Biopsy Site Marker consists of a disposable applicator containing 3 resorbable polyglycolic acid (PGA) pads. The center pad contains a titanium wireform. A polyethylene glycol (PEG) push pad is present in the proximal end.

Product Code	Wireform Material	Wireform Shape	Probe Compatibility
SMTMT11G	Titanium		11G Mammotome™

The wireform is intended for long-term radiographic marking of the biopsy site. The pads are visible via ultrasound for approximately 3 weeks and are essentially resorbed in approximately 12 weeks.

The syringe-like applicator fits within the Mammotome™ probe to access the biopsy cavity. Pads are deployed from the Applicator through the biopsy device into the biopsy cavity.

### INTENDED USE

The SENOMARK<sup>®</sup> Biopsy Site Marker is intended to radiographically and sonographically mark breast tissue during a percutaneous breast biopsy procedure.

### CONTRAINDICATIONS

None known

### WARNINGS

- Patients with a known hypersensitivity to the materials listed in the device description may suffer an allergic reaction to this implant.
- The SENOMARK<sup>®</sup> is not recommended for use in patients with breast implants.
- Do not use in the presence of infection.
- Avoid the use of excessive force during removal of the applicator to prevent breakage of the applicator tip.
- This device has been designed for single use only. Reusing this medical device bears the risk of cross-patient contamination as medical devices – particularly those with long and small lumina, joints, and/or crevices between components – are difficult or impossible to clean once body fluids or tissues with potential pyrogenic or microbial contamination have had contact with the medical device for an indeterminable period of time. The residue of biological material can promote the contamination of the device with pyrogens or microorganisms which may lead to infectious complications.
- Do not resterilize. After resterilization, the sterility of the product is not guaranteed because of an indeterminable degree of potential pyrogenic or microbial contamination which may lead to infectious complications. Cleaning, reprocessing and/or resterilization of the present medical device increases the probability that the device will malfunction due to potential adverse effects on components that are influenced by thermal and/or mechanical changes.

### PRECAUTIONS

- The device should only be used by physicians trained in percutaneous biopsy procedures.
- Do not use this product if the sterile barrier has been previously opened or if the package is damaged.
- This product is provided sterile, and is intended for single use only. Do Not Resterilize.
- Do not kink the flexible tube.
- Maintain correct alignment of the yellow indicator key (C in Figure 1) with the white arrow of the biopsy probe when dispensing pads.
- Ensure that all pads are dispensed.
- Store at temperatures below 25°C (77°F).

### COMPLICATIONS

Potential complications that may be associated with the use of the SENOMARK<sup>®</sup> are similar to those associated with the use of other biopsy marking devices.

### DIRECTIONS FOR USE

Refer to Figure 1 for Applicator diagram.

- A. Plunger    C. Yellow indicator key  
B. Handle    D. Applicator aperture


Figure 1 - Biopsy Site Marker

1. Perform a "dry tap" to ensure that the sample notch is clear of tissue.
2. Using standard aseptic technique, remove the applicator from the package and inspect for damage. Remove tip cover.
3. Insert the applicator into the biopsy device. Advance the applicator to align the yellow indicator key (C) with the white arrow on the thumbwheel of the biopsy probe. Refer to Figure 2.


Figure 2 - Applicator alignment in Probe

4. Fully seat the applicator in the biopsy probe. Verify proper alignment of the applicator aperture with the biopsy chamber by ensuring proper alignment of the yellow indicator key (C) with the white arrow on the thumbwheel of the biopsy probe (Figure 2).
5. Immediately deploy all pads by advancing the plunger (A).  
**Note:** Avoid pulling the applicator handle during deployment which could cause deployment failure.
6. Rotate the biopsy device cannula 180° to position the sample chamber away from the dispensed pads.
7. Remove the applicator from the biopsy device.  
**Caution:** If resistance is felt while removing the applicator, remove the entire probe/applicator assembly from the patient. Failure to do so may result in breakage of the applicator tip.
8. Close the sample chamber and remove the biopsy device, per manufacturer's instructions.
9. Dispose of the applicator properly.
10. Confirm final marker position with imaging.

### WARRANTY

Bard Peripheral Vascular warrants to the first purchaser of this product that this product will be free from defects in materials and workmanship for a period of one year from the date of first purchase and liability under this limited product warranty will be limited to repair or replacement of the defective product, in Bard Peripheral Vascular's sole discretion or refunding your net price paid. Wear and tear from normal use or defects resulting from misuse of this product are not covered by this limited warranty.

**TO THE EXTENT ALLOWABLE BY APPLICABLE LAW, THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL BARD PERIPHERAL VASCULAR BE LIABLE TO YOU FOR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM YOUR HANDLING OR USE OF THIS PRODUCT.**

Some states/countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your state/country.

An issue or revision date and a revision number for these instructions are included for the user's information on the last page of this booklet. In the event 36 months have elapsed between this date and product use, the user should contact Bard Peripheral Vascular, Inc. to see if additional product information is available.

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

**Do Not Use if Temperature Indicator is Black**

Ne pas utiliser si l'indicateur thermique est noir.  
 Bei schwarz gefärbter Temperaturanzeige nicht verwenden  
 Non utilizzare se l'indicatore della temperatura è nero.  
 No usar si el punto de temperatura tiene color negro.  
 Niet gebruiken wanneer de temperatuursindicator zwart is.  
 Não utilize o indicador de temperatura estiver preto  
 Μη χρησιμοποιείτε εάν ο δείκτης θερμοκρασίας είναι μαύρος.  
 Må ikke anvendes, hvis temperaturpletten er sort  
 Använd inte produkten om temperaturindikatorn är svart.  
 Älä käyttää, jos lämpötilan ilmaisin on musta  
 Skäl ikke brukes hvis temperaturindikatoren er svart.  
 Nie używać, jeśli wskaźnik temperatury jest czarny  
 Ne használnia, ha a hőmérséklet-kijelző feketé!  
 Produkt nepoužívejte, pokud je indikátor teploty černý  
 Sıcaklık Göstergesi Siyah Renkliyse Kullanmayın  
 Запрещается использовать, если индикатор температуры имеет черный цвет.  
 如果溫度指示標籤是黑色的，請勿使用

†

**Gauge Size**

Diamètre  
 Gauge-Maß  
 Dimensioni gauge  
 Tamaño de calibre  
 Naaldidike  
 Calibre  
 Μέγεθος διαμετρήματος  
 Gauge-størrelse  
 Gauge-storlek  
 G-koko  
 Målerstørrelse  
 Rozmiar igły  
 Gauge méret  
 Velikost měřky  
 Kalibre Boyutu  
 Калибр  
 直徑尺寸

**Upper Limit of Temperature**

Limite supérieure de température  
 Obere Temperaturgrenze  
 Temperatura massima  
 Límite superior de temperatura  
 Bovengrens temperatuur  
 Limite superior de temperatura  
 Ανώτατο όριο θερμοκρασίας  
 Øvre temperaturgrense  
 Högsta temperatur  
 Lämpötilan yläraja  
 Øvre temperaturgrense  
 Maksymalna dopuszczalna temperatura  
 A hőmérséklet felső határértéke  
 Horní mez teploty  
 Üst Sıcaklık Limiti  
 Верхний предел температуры  
 溫度上限

**Not Made with Natural Rubber Latex**

Fabriqu e sans latex naturel  
 Nicht aus Naturgummilatem hergestellt.  
 Non prodotto con lattice di gomma naturale  
 Este producto no se fabrica con latex de caucho natural  
 Niet gemaakt met latex van natuurlijk rubber  
 N o cont em latex de borracha natural  
 Δεν κατασκευ ζεται απ  φυσικό ελαστικό λάτεξ  
 Ikke fremstillet med naturlig gummilatem  
 Ej tillverkad med naturgummilatem  
 Valmistuksessa ei ole k ytetty luonnonkumilatemksia  
 Dette produktet er ikke produsert med naturlig gummilatemks  
 Produkt nie jest wytwarzany z lateksu kauczuku naturalnego  
 Nem tartalmaz természetes gumilatemet.  
 Nen  vyroben z p rodn ho pry ov ho latexu  
 Dođal Kau uk Lateksten  retimemiştir  
 He содержит натуральный каучуковый латекс  
 製造未採用天然膠乳

**Sterilized Using Irradiation**

St rilis e par irradiation  
 Durch Bestrahlung sterilisiert  
 Sterilizzato mediante irradiazione  
 Esterilizado mediante irradiaci n  
 Gesteriliseerd met behulp van straling  
 Esterilizado por Irradia o  
 Αποστειρωμένο με χρ ση ακτινοβολίας  
 Steriliseret ved brug af bestr ling  
 Steriliserad med str lning  
 Steriloitu s teilytt m ll   
 Steriliseret ved str lning  
 Produkt sterylizowany przez napromienianie  
 Sug rzs ssal steriliz lt  
 Sterilizov no oz fenim  
 Radyoaktif Işınla Sterilize Edilmiştir  
 Стерилизовано радиацией  
 已使用放射線消毒

**Do Not Use If Package Is Damaged**

Ne pas utiliser si l'emballage est endommag e  
 Nicht verwenden, wenn die Verpackung besch digt ist  
 Non utilizzare se la confezione   danneggiata  
 No usar si el envase est  da ado  
 Niet gebruiken wanneer de verpakking beschadigd is  
 N o utilizar se a embalagem estiver danificada  
 Μη χρησιμοποιείτε το πακ on εάν η συσκευασία του  χει υποστεί ζημι   
 M  ikke bruges, hvis emballagen er beskadiget  
 Anv nd inte produkten om f rpackningen skadats  
 Ei saa k ytt t , jos pakkaus on vaurioitunut  
 Sk l ikke brukes hvis pakningen er skadet  
 Nie stosowa , jeřli opakowanie jest uszkodzone  
 Ne haszn lja fel, ha a csomagol s s r lt  
 Nepoužívejte, pokud je obal pořkozen  
 Paket Hasarlıysa Kullanmayınız  
 He использовать, если упаковка повреждена  
 如果包裝已受損，請勿使用

**Breast Tissue Marker**

Marqueur de tissu mammaire  
 Brustgewebemarker  
 Marcatore tessutale mammario  
 Marcador de tejido mamario  
 Mammaweefselmarkering  
 Marcador de tecido mam rio  
 Δείκτης μαστικής ιστού  
 Brystv vsmark r  
 Br sttv vnadsmark r  
 Rintakudosmarkkeri  
 Brystv vmark r  
 Marker tkanki piersiowej  
 Eml sz vet marker  
 Zna kova  prsn  tk n   
 Meme Dokusu İşaretleyicisi  
 Биологический маркер молочной железы  
 乳房組織標記

**Authorized Representative in the European Community**

Repr sentant autoris  au sein de la Communaut  europ enne  
 Berechtigter Vertreter in der Europ ischen Gemeinschaft  
 Rappresentante autorizzato nella Comunit  Europea  
 Representante autorizado en la Comunidad Europea  
 Gemachtigde vertegenwoordiger binnen de Europese Gemeenschap  
 Representante autorizado na Comunidade Europeia  
 Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα  
 Autoriseret repr sentant i Det Europ iske F llesskab  
 Auktoriserad representant i EG  
 Valtuutettu edustaja Euroopan yhteis ss   
 Autoriseret representant i EU  
 Autoryzowany przedstawiciel na terenie Unii Europejskiej  
 Bejegyzett k pviselet az Eur pai K z ss gben  
 Autorizovaný z stupce pro Evropskou unii  
 Avrurpa Topluļuğu Yetkili Temsilcisi  
 Уполномоченный представитель в Европейском сообществе  
 歐洲共同體的授權代表

**Copyright © 2015 C.R. Bard, Inc. All rights reserved.**

Copyright © 2015 C.R. Bard, Inc. Tous droits réservés.  
Copyright © 2015 C.R. Bard, Inc. Alle Rechte vorbehalten.  
Copyright © 2015 C.R. Bard, Inc. Tutti i diritti riservati.  
Copyright © 2015 C.R. Bard, Inc. Todos los derechos reservados.  
Copyright © 2015 C.R. Bard, Inc. Alle rechten voorbehouden.  
Copyright © 2015 C.R. Bard, Inc. Todos os direitos reservados.  
Πνευματικά δικαιώματα © 2015 C.R. Bard, Inc. Με την επιφύλαξη παντός δικαιώματος.  
Copyright © 2015 C.R. Bard, Inc. Alle rettigheder forbeholdt.  
Copyright © 2015 C.R. Bard, Inc. Eftertryck förbjudes.  
Copyright © 2015 C.R. Bard, Inc. Kaikki oikeudet pidätetään.  
Copyright © 2015 C.R. Bard, Inc. Med enerrett.  
Copyright © 2015 C.R. Bard, Inc. Wszelkie prawa zastrzeżone.  
Copyright © 2015 C.R. Bard, Inc. Minden jog fenntartva.  
Copyright © 2015 C.R. Bard, Inc. Všechna práva vyhrazena.  
Telif Hakki © 2015 C.R. Bard, Inc. Tüm hakları saklıdır.  
Авторские права © 2015 Компания С. Р. Bard, Inc. Все права защищены.  
Copyright © 2015 C.R. Bard, Inc. 保留所有權利。  
© 2015 C.R. Bard, Inc.

**Bard and SenoMark are trademarks and/or registered trademarks of C. R. Bard, Inc. or an affiliate. All other trademarks are the property of their respective owners.**

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 **Manufacturer:**  
**SenoRx, Inc.**  
1625 West 3rd Street  
Tempe, AZ 85281  
USA

TEL: 1-480-894-9515  
1-800-321-4254  
FAX: 1-480-966-7062  
1-800-440-5376  
[www.bardbiopsy.com](http://www.bardbiopsy.com)

**CE**  
**0086**

**EC REP**

**Authorized Representative in  
the European Community**  
**Bard Limited**  
Forest House  
Tilgate Forest Business Park  
Brighton Road, Crawley  
West Sussex  
RH11 9BP, UK

**BARB** *Biopsy*

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Bard SenoMark C. R. Bard, Inc.