STARCHMARK®

STARCHMARK® Breast Tissue Marker

For use with ENCOR® Biopsy Probes

ENGLISH

REF STMKEC10GSS
STMKEC7GSS

INSTRUCTIONS FOR USE

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

DESCRIPTION

The STARCHMARK® Breast Tissue Marker is a sterile, disposable Applicator containing 6 resorbable polysaccharide (starch) pellets, 1 polylactic/polyglycolic acid-based copolymer (PLA/PGA) pellet with an embedded radiopaque wireform, and a polyethylene glycol push pellet.

Product Code	Wireform Material	Wireform Shape	Probe Compatibility
STMKEC10GSS	316L stainless steel		10G ENCOR®
STMKEC7GSS	316L stainless steel		7G ENCOR®

The wireform is intended for long-term radiographic marking of the biopsy site. The polysaccharide (starch) pellets absorb body fluids to help in the control and management of bleeding. The pellets are essentially resorbed in approximately 2 weeks.

The syringe-like applicator fits within the ENCOR® probe to access the biopsy cavity. Pellets are deployed from the Applicator through the biopsy device into the biopsy cavity.

INDICATIONS FOR USE

The STARCHMARK® Breast Tissue Marker is intended to radiographically mark breast tissue during a percutaneous breast biopsy procedure.

CONTRAINDICATIONS

Do not use this device in patients with a known hypersensitivity to the materials listed in the device description as they may suffer an allergic reaction to this implant.

WARNINGS

- Use caution when placing near a breast implant to avoid puncture of the implant capsule.
- Do not use in the presence of infection.
- Avoid the use of excessive force during removal of the Applicator to prevent breakage of the Applicator tip.
- This device has been designed for single use only. Reusing this medical device bears the risk of cross-patient contamination as medical devices – particularly those with long and small lumina, joints, and/or crevices between components – are difficult or impossible to clean once body fluids or tissues with potential pyrogenic or microbial contamination have had contact with the medical device for an indeterminable period of time. The residue of biological material can promote the contamination of the device with pyrogens or microorganisms which may lead to infectious complications.
- Do not resterilize. After resterilization, the sterility of the product is not guaranteed because of an indeterminable degree of potential pyrogenic or microbial contamination which may lead to infectious complications. Cleaning, reprocessing and/or resterilization of the present medical device increases the probability that the device will malfunction due to potential adverse effects on components that are influenced by thermal and/or mechanical changes.

PRECAUTIONS

- This product should only be used by a physician who is competently familiar with the indications, contraindications, limitations, typical findings and possible side effects of breast tissue marker placement.

- Do not use this product if the sterile barrier has been previously opened or if the package is damaged.
- After use, this product may be a potential biohazard. Handle and dispose of in accordance with acceptable medical practice and applicable local, state, and federal laws and regulations.
- Do not kink the flexible tube.
- Maintain correct alignment of the yellow indicator key (D in Figure 1) with the red arrow of the biopsy probe when dispensing pellets.
- Ensure that all pellets are dispensed.
- Published studies with comparably sized 316L stainless steel biopsy markers have shown no marker movement and insignificant marker heating when tested in an MRI System with a 1.5T field strength.
- Store at temperatures below 25°C (77°F).
- Do not use if temperature dot is black.
- Although starch (polysaccharide) has known hemostatic properties, the user should continue to employ standard methods for obtaining hemostasis following the biopsy procedure.

POTENTIAL COMPLICATIONS

Complications may occur at any time during or after the procedure. Potential complications of marker placement may include, but are not limited to: hematoma, hemorrhage, infection, adjacent tissue injury and pain.

HOW SUPPLIED

The STARCHMARK® Breast Tissue Marker is provided sterile, and is intended for single use only. Do not reuse. Do not resterilize.

EQUIPMENT REQUIRED

The following equipment is required to place the breast tissue marker:

- Appropriate imaging modality and accessories
- Appropriate biopsy probe
- Surgical gloves and drapes
- Local anesthetic
- Other equipment as necessary

DIRECTIONS FOR USE

Refer to Figure 1 for Applicator diagram.


Figure 1 – Applicator

1. Perform a "dry tap" to ensure that the sample notch is clear of tissue.
2. Using standard aseptic technique, remove the Applicator from the package and inspect for damage. Remove Tip Protector (A).
3. Press the "Marker or Anesthetic" button on the Control Module.
4. Remove the Port Cap from of the Sample Chamber of the probe.

5. Insert Tip Guide (B) into Sample Chamber. Ensure that the Tip Guide is fully inserted. Refer to Figure 2.


Figure 2: Insert Tip Guide into Sample Chamber

6. Remove the Tip Protector from the Applicator.
7. Insert the Applicator through the Tip Guide (B) into the biopsy probe.
Note: A small amount of resistance may be felt as the Applicator Handle (C) passes through the Tip Guide.
8. When the Applicator contacts the proximal end of the probe, rotate the Applicator, if necessary, and continue to advance the handle.
9. Advance the Applicator to align the Yellow Indicator (D) with the red arrow of the biopsy probe. Refer to Figure 3.


Figure 3: Applicator alignment in probe

Caution: To properly deploy pellets, the Applicator Yellow Indicator (D) must be fully seated against the probe and aligned with the red arrow. See Figure 4.


Figure 4: Fully seated Applicator in probe

10. Immediately deploy all pellets by advancing the Plunger (E).
Note: Avoid pulling the Applicator Handle during deployment which could cause deployment failure.
11. Rotate the biopsy device 180° to position the Sample Chamber away from the dispensed pellets.
12. Remove the Applicator from the biopsy device.
Caution: If resistance is felt while removing the Applicator, remove the entire probe/Applicator assembly from the patient. Failure to do so may result in breakage of the Applicator Tip.
13. Dispose of the Applicator properly.
14. Confirm final marker position with imaging.

WARRANTY

Bard Peripheral Vascular warrants to the first purchaser of this product, that this product will be free from defects in materials and workmanship for a period of one year from the date of first purchase and liability under this limited product warranty will be limited, to repair or replacement of the defective product, in Bard Peripheral Vascular's sole discretion, or refunding your net price paid. Wear and tear from normal use or defects resulting from misuse of this product are not covered by this limited warranty.

TO THE EXTENT ALLOWABLE BY APPLICABLE LAW, THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL BARD PERIPHERAL VASCULAR BE LIABLE TO YOU FOR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM YOUR HANDLING OR USE OF THIS PRODUCT.

Some states/countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your state/country.

An issue or revision date and revision number for these instructions are included for the user's information on the last page of this booklet. In the event 36 months have elapsed between this date and product use, the user should contact Bard Peripheral Vascular to see if additional product information is available.

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**Bard and StarchMark are trademarks and/or registered trademarks of C. R. Bard, Inc. or an affiliate. All other trademarks are the property of their respective owners.**
Copyright ©2015 C. R. Bard, Inc. All rights reserved.
Copyright ©2015 C. R. Bard, Inc. Tous droits réservés.
Copyright ©2015 C. R. Bard, Inc. Alle Rechte vorbehalten.
Copyright ©2015 C. R. Bard, Inc. Tutti i diritti riservati.
Copyright ©2015 C. R. Bard, Inc. Todos los derechos reservados.