SENO^{MARK}
SENO^{MARK}® Biopsy Site Marker
For use with 9 Gauge ATEC™

ENGLISH

REF SMTSU9G

Instructions For Use

Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

DESCRIPTION

The Bard® Inc. **SENO^{MARK}**® Biopsy Site Marker consists of a disposable applicator containing:

- 3 resorbable polyglycolic acid (PGA) pads. The center pad contains a titanium wireform.
- 1 polyethylene glycol (PEG) plug at the proximal end.

Product Code	Wireform Material	Wireform Shape	Probe Compatibility
SMTSU9G	Titanium		9G ATEC™

The wireform is intended for long-term radiographic marking of the biopsy site. The pads are visible via ultrasound for approximately 3 weeks and are essentially resorbed in approximately 12 weeks.

The syringe-like applicator fits within the 9 gauge ATEC™ to access the biopsy cavity.¹ Pads are deployed from the applicator through the biopsy device and into the biopsy cavity.

INTENDED USE

The **SENO^{MARK}**® is intended to radiographically and sonographically mark breast tissue during a percutaneous breast biopsy procedure.

CONTRAINDICATIONS

This device is not intended for use except as indicated.

WARNINGS

- **Patients with a known hypersensitivity to the materials listed in the device description may suffer an allergic reaction to this implant.**
- **The **SENO^{MARK}**® is not recommended for use in patients with breast implants.**
- **Do not use in the presence of infection.**
- **Avoid the use of excessive force during removal of the applicator to prevent breakage of the applicator tip.**
- **This device has been designed for single use only. Reusing this medical device bears the risk of cross-patient contamination as medical devices— particularly those with long and small lumina, joints, and/or crevices between components – are difficult or impossible to clean once body fluids or tissues with potential pyrogenic or microbial contamination have had contact with the medical device for an indeterminable period of time. The residue of biological material can promote the contamination of the device with pyrogens or microorganisms which may lead to infectious complications.**
- **Do not resterilize. After resterilization, the sterility of the product is not guaranteed because of an indeterminable degree of potential pyrogenic or microbial contamination which may lead to infectious complications. Cleaning, reprocessing and/or resterilization of the present medical device increases the probability that the device will malfunction due to potential adverse effects on components that are influenced by thermal and/or mechanical changes.**

PRECAUTIONS

- The device should only be used by physicians trained in percutaneous biopsy procedures.
- Do not use this product if the sterile barrier has been previously opened or if the package is damaged.
- This product is provided sterile, and is intended for single use only. **DO NOT RESTERILIZE.**
- Do not kink the flexible tube.
- Maintain correct alignment of the indicators (C & D in Figure 1) with the sample notch when dispensing pads.
- Ensure that all pads are dispensed.
- Store at temperatures below 25° C (77° F).
- Published studied with comparably sized 316L stainless steel biopsy markers have shown no marker movement and insignificant marker heating when tested in an MRI System with a 1.5T field strength.

COMPLICATIONS

Potential complications (e.g. infection) that may be associated with the use of the **SENO^{MARK}**® are the same as those associated with the use of other biopsy marking devices.

HOW SUPPLIED

The **SENO^{MARK}**® is provided sterile, and is intended for single use only.

DIRECTIONS FOR USE (Refer to Figure 1)


Figure 1

1. Make certain that the collection of biopsy specimens has been completed. Ensure that the sample notch has been cleared of all tissue.
2. Inspect the package to ensure that the package integrity has not been compromised. The product is sterile unless the seal is broken.
3. Using standard aseptic technique, remove the **SENO^{MARK}**® device from the package and check for damage. REMOVE TIP PROTECTOR (G).
4. Remove Adaptor (E) if using ATEC™ 0912-20 (12cm L) or 0912-12 (12cm Petite) Biopsy Device. For the ATEC™ 0909-20 (9cm L) or 0909-12 (9cm Petite) Biopsy Device, slide the Adaptor (E) in the direction of the handle until it engages with the stop (D) on the applicator shaft.
5. Following ATEC™ instructions, remove the driver mechanism from the biopsy device cannula, leaving the cannula in the breast.
6. Place the **SENO^{MARK}**® applicator into the ATEC™ biopsy device cannula and advance until the ATEC™ stop (D) or the Adaptor (E) is seated in the cannula.
7. Align the **SENO^{MARK}**® indicators (C & D) with the sample notch in the biopsy device cannula.
8. Immediately deploy the **SENO^{MARK}**® pads by advancing the plunger (A) while maintaining forward pressure on the handle (B), and correct alignment of the applicator with the sample notch. The pads will be fully dispensed when the plunger contacts the handle.
9. Rotate the biopsy device cannula 180 degrees to position the sample notch away from the dispensed pads.
10. Remove the **SENO^{MARK}**® applicator from the biopsy device cannula. Avoid the use of force. (See Caution below.)
11. Close the sample notch and remove the biopsy device, per manufacturer's instructions.
12. Dispose of the applicator properly.
13. Confirm final marker position with imaging.

Caution: If any resistance is encountered while removing the applicator, leave the **SENO^{MARK}**® applicator within the probe and remove the entire probe/applicator assembly. Failure to do so may result in breakage of the applicator tip.

WARRANTY

Bard Peripheral Vascular warrants to the first purchaser of this product that this product will be free from defects in materials and workmanship for a period of one year from the date of first purchase and liability under this limited product warranty will be limited to repair or replacement of the defective product, in Bard Peripheral Vascular's sole discretion or refunding your net price paid. Wear and tear from normal use or defects resulting from misuse of this product are not covered by this limited warranty.

TO THE EXTENT ALLOWABLE BY APPLICABLE LAW, THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL BARD PERIPHERAL VASCULAR BE LIABLE TO YOU FOR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM YOUR HANDLING OR USE OF THIS PRODUCT.

Some states/countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your state/country.

An issue or revision date and a revision number for these instructions are included for the user's information on the last page of this booklet. In the event 36 months have elapsed between this date and product use, the user should contact Bard Peripheral Vascular, Inc. to see if additional product information is available.

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**Do Not Use if Temperature Indicator is Black**

Ne pas utiliser si l'indicateur thermique est noir.
 Bei schwarz gefärbter Temperaturanzeige nicht verwenden
 Non utilizzare se l'indicatore della temperatura è nero.
 No usar si el punto de temperatura tiene color negro.
 Niet gebruiken wanneer de temperatuursindicator zwart is.
 Não utilize o indicador de temperatura estiver preto
 Μη χρησιμοποιείτε εάν ο δείκτης θερμοκρασίας είναι μαύρος.
 Må ikke anvendes, hvis temperaturpletten er sort
 Använd inte produkten om temperaturindikatorn är svart.
 Älä käytä, jos lämpötilan ilmaisin on musta
 Skäl ikke brukes hvis temperaturindikatoren er svart.
 Nie używać, jeśli wskaźnik temperatury jest czarny
 Ne használja, ha a hőmérséklet-kijelző feketé!
 Produkt nepoužívejte, pokud je indikátor teploty černý
 Sıcaklık Göstergesi Siyah Renkliyse Kullanmayın
 Запрещается использовать, если индикатор температуры имеет черный цвет.
 如果溫度指示標籤是黑色的，請勿使用

†

**Gauge Size**

Diamètre
 Gauge-Maß
 Dimensioini gauge
 Tamaño de calibre
 Naaldidikte
 Calibre
 Μέγεθος διαμετρήματος
 Gauge-størrelse
 Gauge-storlek
 G-koko
 Målerstørrelse
 Rozmiar igły
 Gauge méret
 Velikost měřky
 Kalibre Boyutu
 Калибр
 直徑尺寸

**Upper Limit of Temperature**

Limite supérieure de température
 Obere Temperaturgrenze
 Temperatura massima
 Límite superior de temperatura
 Bovengrens temperatuur
 Limite superior de temperatura
 Ανώτατο όριο θερμοκρασίας
 Øvre temperaturgrense
 Högsta temperatur
 Lämpötilan yläraja
 Øvre temperaturgrense
 Maksymalna dopuszczalna temperatura
 A hőmérséklet felső határértéke
 Horní mez teploty
 Üst Sıcaklık Limiti
 Верхний предел температуры
 溫度上限

**Not Made with Natural Rubber Latex**

Fabriqu  sans latex naturel
 Nicht aus Naturgummil t hergestell.
 Non prodotto con lattice di gomma naturale
 Este producto no se fabrica con l t de caucho natural
 Niet gemaakt met latex van natuurlijk rubber
 N o cont m l t de borracha natural
 Δεν κατασκευ ζεται απ  φυσικό ελαστικό λ τεξ
 Ikke fremstillet med naturlig gummil t
 Ej tillverkad med naturgummil t
 Valmistuksessa ei ole k ytetty luonnonkumil teksia
 Dette produktet er ikke produsert med naturlig gummil t
 Produkt nie jest wytwarzany z lateksu kauczuku naturalnego
 Nem tartalmaz természetes gumil t.
 Nen  vyroben z p rodn ho pry ov ho latexu
 Dođal Kau uk Lateksten  retilimemiştir
 He содержит натуральный каучуковый л текс
 製造未採用天然膠乳

**Sterilized Using Irradiation**

St rilis  par irradiation
 Durch Bestrahlung sterilisiert
 Sterilizzato mediante irradiazione
 Esterilizado mediante irradiaci n
 Gesteriliseerd met behulp van straling
 Esterilizado por Irradia o
 Αποστειρωμένο με χρ ση ακτινοβολίας
 Steriliseret ved brug af bestr ling
 Steriliserad med str lning
 Steriloitu s teilytt m ll 
 Steriliseret ved str ling
 Produkt sterylizowany przez napromienianie
 Sug rzs ssal steriliz lt
 Sterilizov no oz a enim
 Radyoaktif Işınıla Sterilize Edilmiştir
 Стерилизовано радиацией
 已使用放射線消毒

**Do Not Use If Package Is Damaged**

Ne pas utiliser si l'emballage est endommag 
 Nicht verwenden, wenn die Verpackung besch digt ist
 Non utilizzare se la confezione   danneggiata
 No usar si el envase est  da ado
 Niet gebruiken wanneer de verpakking beschadigd is
 N o utilizar se a embalagem estiver danificada
 Μη χρησιμοποιείτε το πακ on εάν η συσκευασία του  χει υποστεί ζημι 
 M  ikke bruges, hvis emballagen er beskadiget
 Anv nd inte produkten om f rpackningen skadats
 Ei saa k ytt , jos pakkaus on vaurioitunut
 Sk l ikke brukes hvis pakningen er skadet
 Nie stosowa , jeřli opakowanie jest uszkodzone
 Ne haszn lja fel, ha a csomagol s s r lt
 Nepoužívejte, pokud je obal pořkozen
 Paket Hasarlıysa Kullanmayınız
 He использовать, если упаковка повреждена
 如果包裝已受損，請勿使用

**Breast Tissue Marker**

Marqueur de tissu mammaire
 Brustgewebemarker
 Marcatore tessutale mammario
 Marcador de tejido mamario
 Mammaweefselmarkering
 Marcador de tecido mam rio
 Δείκτης μαστικής ιστού
 Brystv vsmark r
 Br stv vnadsmark r
 Rintakudosmarkkeri
 Brystv vmark r
 Marker tkanki piersiowej
 Eml sz vet marker
 Zna kova  prsn  tk n 
 Meme Dokusu İşaretleyicisi
 Биологический маркер молочной железы
 乳房組織標記

**Authorized Representative in the European Community**

Repr sentant autoris  au sein de la Communaut  europ enne
 Berechtigter Vertreter in der Europ ischen Gemeinschaft
 Rappresentante autorizzato nella Comunit  Europea
 Representante autorizado en la Comunidad Europea
 Gemachtigde vertegenwoordiger binnen de Europese Gemeenschap
 Representante autorizado na Comunidade Europeia
 Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα
 Autoriseret repr sentant i Det Europ iske F llesskab
 Auktoriserad representant i EG
 Valtuutettu edustaja Euroopan yhteis ss 
 Autoriseret representant i EU
 Autoryzowany przedstawiciel na terenie Unii Europejskiej
 Bejegyzett k pviselet az Eur pai K z ss gben
 Autorizovaný z stupce pro Evropskou unii
 Avrurpa Topluļuğu Yetkili Temsilcisi
 Уполномоченный представитель в Европейском сообществе
 歐洲共同體的授權代表

Copyright © 2015 C.R. Bard, Inc. All rights reserved.

Copyright © 2015 C.R. Bard, Inc. Tous droits réservés.
Copyright © 2015 C.R. Bard, Inc. Alle Rechte vorbehalten.
Copyright © 2015 C.R. Bard, Inc. Tutti i diritti riservati.
Copyright © 2015 C.R. Bard, Inc. Todos los derechos reservados.
Copyright © 2015 C.R. Bard, Inc. Alle rechten voorbehouden.
Copyright © 2015 C.R. Bard, Inc. Todos os direitos reservados.
Πνευματικά δικαιώματα © 2015 C.R. Bard, Inc. Με την επιφύλαξη παντός δικαιώματος.
Copyright © 2015 C.R. Bard, Inc. Alle rettigheder forbeholdt.
Copyright © 2015 C.R. Bard, Inc. Eftertryck förbjudes.
Copyright © 2015 C.R. Bard, Inc. Kaikki oikeudet pidätetään.
Copyright © 2015 C.R. Bard, Inc. Med enerrett.
Copyright © 2015 C.R. Bard, Inc. Wszelkie prawa zastrzeżone.
Copyright © 2015 C.R. Bard, Inc. Minden jog fenntartva.
Copyright © 2015 C.R. Bard, Inc. Všechna práva vyhrazena.
Telif Hakkı © 2015 C.R. Bard, Inc. Tüm hakları saklıdır.
Авторские права © 2015 Компания С. Р. Bard, Inc. Все права защищены.
Copyright © 2015 C.R. Bard, Inc. 保留所有權利。
© 2015 C.R. Bard, Inc.

Bard and SenoMark are trademarks and/or registered trademarks of C. R. Bard, Inc. or an affiliate. All other trademarks are the property of their respective owners.

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 **Manufacturer:**
SenoRx, Inc.
1625 West 3rd Street
Tempe, AZ 85281
USA

TEL: 1-480-894-9515
1-800-321-4254
FAX: 1-480-966-7062
1-800-440-5376
www.bardbiopsy.com

CE
0086

EC REP

**Authorized Representative in
the European Community**
Bard Limited
Forest House
Tilgate Forest Business Park
Brighton Road, Crawley
West Sussex
RH11 9BP, UK

BARB *Biopsy*

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Bard SenoMark C. R. Bard, Inc.