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## Instructions for Use:

**Rx Only: Federal (USA) law restricts this device to sale by or on the order of a physician.**

### A. General Information and Device Description:

The Bard\* Navarre\* Universal Drainage Catheter with Nitinol includes a drainage catheter and an introduction set comprised of a cannula and a needle trocar. The catheter is 30 cm long and is available in diameters from 6 to 12 French. The 6 and 8 French diameter catheters are also available in a 15 cm long configuration. Three tip configurations are available: J-curve, pigtail, and locking pigtail. The catheter body is embedded with Nitinol reinforcement.

The Bard\* Navarre\* OPTI-DRAIN\* Multi-Use Drainage Catheter includes a drainage catheter and an introduction set comprised of a cannula and a needle trocar. The catheter is 30 cm long and is available in diameters from 8 to 12 French. Three tip configurations are available: J-curve, pigtail, and locking pigtail.

The recommended maximum indwelling time for Navarre Drainage Catheters is 10 days. This recommendation is based on results of preclinical in vivo implant testing where the maximum implant duration was 10 days. Bard makes no claims for this product regarding biocompatibility beyond 10 days.

### B. Indications for Use:

The Bard\* Navarre\* Universal Drainage Catheter with Nitinol and the Bard\* Navarre\* OPTI-DRAIN\* Multi-Use Drainage Catheter are intended to be used for abscess, cyst, and other general purpose drainage applications.

### C. Contraindications:

This device is not intended for intravascular use. Drainage catheters should not be used in the absence of a safe route to fluid collection. Coagulation deficiencies should also be addressed prior to drainage catheter use.

### D. Warnings:

- 1. Before using, inspect the drainage catheter for damage. If the catheter is damaged, DO NOT USE.**
- 2. After use, the product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practices and applicable local, state and federal laws and regulations.**
- 3. The device has not been tested or shown to be safe for subjects undergoing MRI.**

### E. Precautions:

- 1. For suspected or confirmed sepsis, ensure that adequate antibiotic coverage is initiated prior to use.**
- 2. Intended for single use. DO NOT REUSE. Reuse and/or repackaging may create a risk of patient or user infection, compromise the structural integrity and/or essential material and design characteristics of the device, which may lead to device failure, and/or lead to injury, illness or death of the patient.**
- 3. Not recommended for long term use. Exchange as needed to maintain patency and minimize risk of leakage.**
- 4. Do not use past expiry date.**
- 5. Do not use excessive force when inserting cannula to straighten catheter. Catheter tip may be straightened manually to aid in insertion of cannula.**
- 6. If trocar is to be used, ensure that hang tag between trocar and cannula hubs remains in place during insertion.**

### F. Potential Complications:

Potential complications of percutaneous drainage catheter procedures are site specific and may consist of hematoma, hemorrhage, infection, adjacent tissue injury, and pain.

#### G. Directions for Use:

1. Select a puncture site. Use diagnostic imaging to locate a safe catheter tract.
2. Anesthetize the puncture site as needed.
3. Prepare the drainage catheter for insertion by advancing the cannula into the drainage catheter. Lock the cannula hub into position. **CAUTION: Do not use excessive force when inserting cannula to straighten catheter. Catheter tip may be straightened manually to aid in insertion of cannula.**
4. For Seldinger technique, remove the stylet from the cannula and advance the catheter over the guidewire using diagnostic imaging to direct the catheter into fluid collection site. Remove the cannula and guidewire when the catheter is in proper position.  
**CAUTION: If trocar is to be used, ensure that hang tag between trocar and cannula hubs remains in place during insertion.**
5. For trocar technique, advance the stylet into the cannula and lock the hub in position. Make a small incision slightly larger than the diameter of the drainage catheter. Use diagnostic imaging to direct the catheter into fluid collection site. Remove the cannula and trocar stylet when the catheter is in proper position.
6. Aspirate the fluid collection and confirm position with diagnostic imaging.
7. Secure the catheter to the skin.
8. Using aseptic technique, securely attach accessories, as needed, and inspect for adequate drainage. Use gravity or low suction to provide continued drainage.

#### H. Suggested Method to Remove or Replace the Catheter:

1. Disconnect any connections from drainage catheter.
2. Pull the catheter out gently.
3. If access to the drainage site is to be maintained, remove any connections if applicable, as described in methods #1 through #2. Advance a guidewire past the distal tip of the catheter. Withdraw catheter over guidewire.

An issued or revision date and revision number for these instructions are included for the users information below. In the event two years have elapsed between this date and product use, the user should contact Bard Access Systems, Inc. to see if additional product information is available.

#### Revised Date: September 2010

\*Bard, Navarre, and Opti-Drain are trademarks and/or registered trademarks of C. R. Bard, Inc.

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© 2010 C. R. Bard, Inc. All rights reserved.
Bard Limited