ENCOR® Directional Vacuum-Assisted Biopsy Device and Driver

ENGLISH

INSTRUCTIONS FOR USE

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

DESCRIPTION

The BARD ENCOR® Breast Biopsy Device (probe) is indicated for use in acquiring tissue for diagnosing breast abnormalities. The ENCOR® Breast Biopsy Probe and ENCOR® Breast Biopsy Driver are sold and packaged separately. The ENCOR® Breast Biopsy Probe and Breast Biopsy Driver are used together with one of the ENCOR® Breast Biopsy Systems which include ENCOR®, ENCOR® ENSPIRE, and ENCOR® ULTRA. These systems provide vacuum and system control. Refer to the Instructions for Use (IFU) for the specific ENCOR® Breast Biopsy System being used.

The ENCOR® Breast Biopsy Probe contains a sharp trocar tip at the distal end for insertion and oscillating tubular cutter for tissue acquisition. A detachable tissue collection chamber is located on the proximal end of the device. The tissue collection chamber incorporates both a suction line for acquisition and transport of tissue samples and a saline rinse line for optional rinse of acquired samples in the chamber. A syringe adapter is provided for introducing anesthetic to the biopsy site. BARD® Breast Tissue Markers which are sold and packaged separately may be inserted through the proximal end of the device.

The ENCOR® Breast Biopsy Driver is designed as a handheld unit for ultrasound-guided breast biopsies and for mounting on a stereotactic platform using BARD adapters. These adapters are sold and packaged separately.

INTENDED USE

The ENCOR® Breast Biopsy Probe and ENCOR® Breast Biopsy Driver are indicated to acquire tissue for diagnostic sampling of breast abnormalities. The instrument is intended to provide breast tissue for histologic examination with partial or complete removal of the imaged abnormality.

The extent of histologic abnormality cannot be reliably determined from its mammographic appearance. Therefore the extent of removal of the imaged evidence of an abnormality does not predict the extent of removal of a histologic abnormality, e.g., malignancy. When the sampled abnormality is not histologically benign, it is essential that the tissue margins be examined for completeness of removal using standard surgical procedures.

CONTRAINDICATIONS

This device is not intended for use except as indicated. In addition, it is not intended for use under conditions where breast tissue excision is contraindicated.

WARNINGS

- Care must be taken when positioning the device trocar tip near the chest wall or skin margin.
- This device has been designed for single use only. Reusing this medical device bears the risk of cross-patient contamination as medical devices – particularly those with long and small lumina, joints, and/or crevices between components – are difficult or impossible to clean once body fluids or tissues with potential pyrogenic or microbial contamination have had contact with the medical device for an indeterminable period of time. The residue of biological material can promote the contamination of the device with pyrogens or microorganisms which may lead to infectious complications.
- Do not resterilize. After resterilization, the sterility of the product is not guaranteed because of an indeterminable degree of potential pyrogenic or microbial contamination which may lead to infectious complications. Cleaning, reprocessing and/or resterilization of the present medical device increases the probability that the device will malfunction due to potential adverse effects on components that are influenced by thermal and/or mechanical changes.

- Di(2-ethylhexyl)phthalate (DEHP) is a plasticizer used in some polyvinyl chloride medical devices. DEHP has been shown to produce a range of adverse effects in experimental animals, notably liver toxicity and testicular atrophy. Although the toxic and carcinogenic effects of DEHP have been well established in experimental animals, the ability of this compound to produce adverse effects in humans is controversial. There is no evidence that neonates, infants, pregnant and breast feeding women exposed to DEHP experience any related adverse effects. However, a lack of evidence of causation between DEHP-PVC and any disease or adverse effect does not mean that there are no risks.

PRECAUTIONS

- A scalpel should be used to incise the skin prior to inserting the trocar tip.
- The ENCOR® Probe is provided sterile and is intended for single use only. Do not resterilize.
- The ENCOR® Driver is reusable and is supplied non-sterile. Do not autoclave or immerse in liquid.

COMPLICATIONS

Complications that may be associated with the use of the ENCOR® Biopsy Device and Driver are the same as those associated with the use of other biopsy devices. These may include injury to the skin, blood vessels, muscles, organs, bleeding, hematoma and infection.

HOW SUPPLIED

The ENCOR® Biopsy Probe is supplied sterile and is intended for single use only.

The ENCOR® Driver is reusable and is supplied non-sterile.

DIRECTIONS FOR USE

1. Prepare the ENCOR® Breast Biopsy System (ENCOR®, ENCOR® ENSPIRE, ENCOR® ULTRA) according to the respective Instructions for Use (IFU).
2. Turn on the ENCOR® Breast Biopsy System (ENCOR®, ENCOR® ENSPIRE, ENCOR® ULTRA) and allow it to initialize. The ENCOR® Breast Biopsy System will indicate when the ENCOR® Breast Biopsy Driver initialization was successful and that the system is ready for the ENCOR® Breast Biopsy Probe installation. Note: the ENCOR® Breast Biopsy Driver will not initialize with the ENCOR® Breast Biopsy Probe installed.
3. For stereotactic table use, refer to the instructions for use provided with the appropriate table adapter for proper installation of the ENCOR® Breast Biopsy Driver.
4. Inspect the ENCOR® Breast Biopsy Probe package to ensure that the package integrity has not been compromised. The product is sterile unless the package has been previously opened or damaged. Do not use the product if the sterility has been compromised.
5. Using standard aseptic technique, remove the ENCOR® Breast Biopsy Probe from the package and check for damage. Do not use the probe if damaged. Install the ENCOR® Breast Biopsy Probe to the ENCOR® Breast Biopsy Driver by sliding the distal end into place and then pressing down on the proximal end to lock (see Figure 1). Connect the vacuum tubing and optional saline rinse connection to the ENCOR® Breast Biopsy System (ENCOR®, ENCOR® ENSPIRE, ENCOR® ULTRA) according to the appropriate IFU. For stereotactic use, remove the plastic cover from the tip.

Note: To install the optional tissue collection tray cover (ECTT), depress the finger tabs to remove the tissue collection tray from the tissue collection chamber of the ENCOR® Breast Biopsy Probe. Insert the tissue collection tray cover onto the top of the tissue collection tray. Replace the tissue collection tray into the ENCOR® Breast Biopsy Probe. (See Figure 2)

Figure 1: Attach ENCOR® Breast Biopsy Probe to ENCOR® Breast Biopsy Driver

6. Calibrate the ENCOR® Breast Biopsy Probe for handheld or stereotactic use according to the ENCOR® Breast Biopsy System Instructions for Use. Exercising caution, remove the tip protector from the distal end of the ENCOR® Breast Biopsy Probe.
7. Select desired sampling options according to the ENCOR® Breast Biopsy System (ENCOR®, ENCOR® ENSPIRE, ENCOR® ULTRA) Instructions for Use.
8. Use standard techniques to anesthetize the area.
9. Using a scalpel (not provided), make a skin incision and position the tip of the ENCOR® Breast Biopsy Probe subcutaneously.
10. Advance the ENCOR® Breast Biopsy Device and position at the target site.
11. To initiate the sampling sequence, utilize either the "Sample" button on the ENCOR® Breast Biopsy Driver or "Sample" switch on the ENCOR® Breast Biopsy Foot Pedal as appropriate. While firmly holding the ENCOR® Breast Biopsy Driver to maintain the aperture at the target site, press the "Sample" button or foot switch again to initiate the automated tissue acquisition cycle. A tissue sample will be automatically transported to the tissue collection chamber. If the rinse sample is enabled, each sample will be rinsed with saline. Repeat this step as needed to obtain sufficient tissue samples.

To evacuate fluid or blood, press either the VAC button on the ENCOR® Breast Biopsy Driver or VAC button on the foot pedal.

Note: If the ENCOR® Breast Biopsy System (ENCOR®, ENCOR® ENSPIRE, ENCOR® ULTRA) alarm tone sounds during the procedure, refer to the ENCOR® Breast Biopsy System display screen for information on how to proceed.

12. See IFU of appropriate ENCOR® Breast Biopsy System (ENCOR®, ENCOR® ENSPIRE, ENCOR® ULTRA) for instructions on utilizing optional anesthetic delivery and sample saline flush.
13. To remove tissue samples, press both finger tabs and separate the rear tissue collection tray from the tissue collection chamber. If a tissue tray cover is used, remove tissue tray cover to retrieve samples from tissue collection tray. The tissue collection tray may be repositioned to acquire additional samples from the breast.
14. At the conclusion of the procedure, remove the ENCOR® Breast Biopsy Probe from the breast. Disconnect ENCOR® Breast Biopsy Probe from ENCOR® Breast Biopsy Driver and clean up ENCOR® Breast Biopsy System (ENCOR®, ENCOR® ENSPIRE, ENCOR® ULTRA) according to appropriate IFU.

Note: Clean the ENCOR® Breast Biopsy Driver by wiping all external surfaces with a cloth soaked with Dispatch Hospital cleaner or other hospital approved cleaner/disinfectant.

Note: The ENCOR® Breast Biopsy Probe and components are for single use only. After use, these products may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL BARD PERIPHERAL VASCULAR BE LIABLE TO YOU FOR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM YOUR HANDLING OR USE OF THIS PRODUCT.

Some states/countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your state/country.

An issue or revision date and a revision number for these instructions are included for the user's information on the last page of this booklet. In the event 36 months have elapsed between this date and product use, the user should contact Bard Peripheral Vascular, Inc. to see if additional product information is available.

Assembled in Thailand.


Figure 2: ENCOR® Breast Biopsy Probe

WARRANTY

Bard Peripheral Vascular warrants to the first purchaser of this product that this product will be free from defects in materials and workmanship for a period of one year from the date of first purchase and liability under this limited product warranty will be limited to repair or replacement of the defective product, in Bard Peripheral Vascular's sole discretion or refunding your net price paid. Wear and tear from normal use or defects resulting from misuse of this product are not covered by this limited warranty.

TO THE EXTENT ALLOWABLE BY APPLICABLE LAW, THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO,

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**Contains or presence of Phthalates**

Contient ou présence de phtalates
Contiene o sono presenti ftalati
Contiene ftalatos o se ha determinado la presencia de ftalatos
Enthält Phthalate
Bevat ftalaten of ftalaten aanwezig
Contém ou presença de ftalatos
Περιέχει φθαλικές ενώσεις ή ίχνη αυτών
Indeholder eller tilstedeværelse af phthalater
Innehåller ftalater
Sisältää ftalaatteja tai ftalaatteja on läsnä
Inneholder ftalater eller ftalater er til stede
Zawartość lub obecność ftalatów
Ftalátot tartalmaz
Obsahuje ftaláty nebo jejich stopy
Ftalatlar içerir veya vardır
Содержание или присутствие фталатов
包含或存在鄰苯二甲酸鹽
프탈레이트류가 포함되어 있거나 검출됨

**Sterilized Using Irradiation**

Stérilisé par irradiation
Durch Bestrahlung sterilisiert
Sterilizzato mediante irradiazione
Esterilizado mediante irradiación
Gesteriliseerd met behulp van straling
Esterilizado por Irradiação
Αποστειρωμένο με χρήση ακτινοβολίας
Steriliseret ved brug af bestråling
Steriliserad med strålning
Steriloitu säteilyttämällä
Sterilisert ved stråling
Produkt sterylizowany przez napromienianie
Sugárzással sterilizált
Sterilizováno ozářením
Radyoaktif Işınla Sterilize Edilmiştir
Стерилизовано радиацией
已使用放射線消毒
방사선을 사용한 멸균 처리

**Do Not Use if the Product Sterile Barrier System or its Packaging is Compromised**

Ne pas utiliser si le système de barrière stérile du produit ou son emballage sont ou paraissent endommagés
Bei beeinträchtigtem Sterilschutzsystem oder beschädigter Verpackung das Produkt nicht verwenden
Non utilizzare se il sistema della barriera sterile del prodotto o la confezione sono compromessi
No utilizar si el sistema de barrera estéril del producto o su envase están dañados
Niet gebruiken wanneer het steriele-barrièresysteem van het product of de verpakking is aangetast
Não utilizar se o sistema de barreira estéril do produto ou respectiva embalagem estiverem comprometidos
Μη χρησιμοποιείτε το προϊόν εάν έχει παραβιαστεί το σύστημα φραγμού αποστείρωσης ή η συσκευασία του
Må ikke anvendes, hvis produktets sterile barrieresystem eller emballagen er kompromitteret
Använd inte om produktens sterila barriär eller förpackning är skadad
Ei saa käyttää, jos tuotteen steriilisuojaus tai pakkaus on vaurioitunut
Skal ikke brukes hvis produktets sterile barrieresystem eller pakningen er forringet
Nie używać, jeśli naruszono sterylność produktu lub jego opakowanie
Ne használja, ha a termékét védő steril lezárási rendszer vagy a csomagolás sérült!
Pokud je narušena sterilní bariéra nebo poškozen obal, výrobek nepoužívejte
Ürünün Steril Bariyer Sistemi ya da Ambalajı Bozulmuşsa Ürünü Kullanmayın
Запрещается применять изделие, если его стерильная упаковка или внешняя упаковка повреждены
如果產品消毒屏障系統或包裝受損，請勿使用
제품의 멸균 장벽 시스템 또는 포장에 손상된 경우 사용하지 마십시오.

**Authorised Representative in the European Community**

Représentant autorisé au sein de la Communauté européenne
Bevollmächtigter in der Europäischen Gemeinschaft
Rappresentante autorizzato nella Comunità Europea
Representante autorizado en la Comunidad Europea
Gemachtigde vertegenwoordiger binnen de Europese Gemeenschap
Representante autorizado na Comunidade Europeia
Εξουσιοδοτημένος Αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα
Autoriseret repræsentant i Det Europæiske Fællesskab
Auktoriserad representant i EG
Valtuutettu edustaja Euroopan yhteisössä
Autorisert representant i EU
Autoryzowany przedstawiciel w Unii Europejskiej
Bejegyzett képviselő az Európai Közösségben
Autorizovaný zástupce pro Evropské společenství
Avrupa Topluluğu Yetkili Temsilcisi
Уполномоченный представитель в Европейском сообществе
歐洲共同體的授權代表
유럽 공동체 내의 인증받은 대리업체

Do not dispose of the EnCor System by placing into trash receptacles

Ne pas jeter de l'EnCor système par la mise en récipients à ordures
Das EnCor-System nicht in Müllbehälter entsorgen
No deseché el sistema EnCor colocándolo en recipientes de basura
Non smaltire il sistema EnCor nella pattumiera
Het EnCor systeem dient niet in vuilniscontainer te worden verwerkt
Não elimine o sistema EnCor em recipientes para lixo doméstico
Μην απορρίπτετε το σύστημα EnCor σε κάδους απορριμάτων
EnCor-systemet må ikke bortskaffes ved at lægge det i en skraldspannd
EnCor-systemet får inte bortskaffas i vanliga sopbehållare
Älä hävittä EnCor-järjestelmää panemalla se roska-astiaan
Ikke avhend EnCor-systemet ved å legge det i søppelkasser
Nie należy wyrzucać systemu EnCor przez umieszczenie go w pojemnikach na śmieci
Az EnCor rendszer hulladékkezelésekor ne dobja azt a szemetes tartályba
Neodkládejte systém EnCor do běžného odpadu
EnCor Sistemini çöp kovalarına koyarak atmayın
Не выбрасывайте систему EnCor в мусорные контейнеры
棄置 EnCor 系統時，請勿放在垃圾桶
EnCor 시스템을 쓰레기통에 버리지 마십시오

Copyright © 2015 C. R. Bard, Inc. All rights reserved.

Copyright © 2015 C. R. Bard, Inc. Tous droits réservés.
Copyright © 2015 C. R. Bard, Inc. Alle Rechte vorbehalten.
Copyright © 2015 C. R. Bard, Inc. Tutti i diritti riservati.
Copyright © 2015, C. R. Bard, Inc. Reservados todos los derechos.
Copyright © 2015 C. R. Bard, Inc. Alle rechten voorbehouden.
Copyright © 2015, C. R. Bard, Inc. Todos os direitos reservados.
Πνευματικά δικαιώματα © 2015 C. R. Bard, Inc. Με την επιφύλαξη παντός δικαιώματος.
Copyright © 2015 C. R. Bard, Inc. Alle rettigheder forbeholdes.
Copyright © 2015 C. R. Bard, Inc. Alla rättigheter förbehållna.
Copyright © 2015 C. R. Bard, Inc. Kaikki oikeudet pidätetään.
Copyright © 2015 C. R. Bard, Inc. Alle rettigheter forbeholdt.
Copyright © 2015 C. R. Bard, Inc. Wszelkie prawa zastrzeżone.
Copyright © 2015, C. R. Bard, Inc. Minden jog fenntartva.
Copyright © 2015 C. R. Bard, Inc. Všechna práva vyhrazena.
Yapıt hakkı © 2015 C. R. Bard, Inc. Tüm hakları saklıdır.
Авторские права © 2015 г., C. R. Bard, Inc. Все права защищены.
Copyright © 2015 C. R. Bard, Inc. 保留所有權利。
저작권 © 2015 C. R. Bard, Inc. 모든 권리 보유.

Bard, EnCor, EnCor Enspire, EnCor Ultra, Enspire and Ultra are trademarks and/or registered trademarks of C.R. Bard, Inc. or an affiliate

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CE
0086

EC REP

**Authorized Representative in
the European Community**

Bard Limited
Forest House
Tilgate Forest Business Park
Brighton Road, Crawley
West Sussex
RH11 9BP, UK

BARDD *Biopsy*

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SenoRx, Inc.